V 300 S / VS 300 S



EN Installation and operating instructions





Contents

				As	sem	oly	
				6	Requ	irements	22
lm	porta	nt information			6.1	Installation/setup room	22
1	-	this document	3		6.2	Setup options	22
'	1.1	Warnings and symbols	3		6.3	Pipe materials	22
	1.2	Copyright information	4		6.4	Hose materials	22
					6.5	Information about electrical con-	
2		/	4			nections	22
	2.1	Intended purpose	4		6.6	Information about connecting	
	2.2	Intended use	4			cables	22
	2.3	Improper use	4	7	Syste	em components	23
	2.4	Systems, connection with other	_		7.1	OroCup	23
		devices	5		7.2	Condensate separator	23
	2.5	General safety information	5		7.3	Rinsing unit	23
	2.6	Specialist personnel	5		7.4	Flow accelerator	23
	2.7	Notification requirement of seri-	E		7.5	Amalgam separator	23
	0.0	ous incidents	5 5		7.6	Surge tank	24
	2.8	Electrical safety	6		7.7	Bacteria filter	24
	2.9	Only use original parts			7.8	Noise reduction	24
	2.10	Transport	6 6	8	Instal	llation	25
	2.11	Disposai	O		8.1	Mounting the condensation sep-	
E	L					arator on the V 300 S	25
Œ					8.2	Installation and routeing of hoses	
Pr	nduct	description				and pipes	25
		•	7		8.3	Rinsing unit water connections.	28
3	3.1	iew	7 8	9	Elect	rical connections	29
	3.1	Optional items	8	10	Comr	missioning	29
	3.3	Consumables	8			3	
	3.4	Wear parts and replacement	O				
	0.4	parts	8				
4	Toohn	ical data	9	Us	age		
4	4.1	V 300 S	9		-	fection and cleaning	30
	4.2	VS 300 S	13	• • •	11.1	Suctioning water	31
	4.3	Type plate	17		11.2	Disinfection	31
	4.4	Evaluation of conformity	17		11.3	Cleaning	32
_		•		10		-	34
5		ion	18	12	12.1	tenance	34
	5.1	V 300 S	18		12.1		
	5.2	VS 300 S	20		12.2	VS 300 S	34
					12.3	parts	34
					12.4	Cleaning the protective strainer.	35
					12.4	Intensive cleaning	35
					12.0	intonsive clearing	00



Troubleshooting

13	Tips for operators and service technicians	37
14	Transporting the unit	39



Appendix

15	Handover record	40
16	Country representatives	41

Important information

About this document

These installation and operating instructions represent part of the unit.



The manufacturer and the distributor will not offer any guarantee or accept any liability for the safe operation and the safe functioning of the unit if the instructions and information in these installation and operating instructions are not complied with.

The German version of the installation and operating instructions is the original manual. All other languages are translations of the original manual. These installation and operating instructions apply to:

V 300 S

REF: 7119-01: 7119-01/002: 7119-02: 7119-02/002

VS 300 S

REF: 7122-01; 7122-01/002; 7122-01/021; 7122-02; 7122-02/002; 7122-04; 7122-04/002

1.1 Warnings and symbols

Warnings

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The following warning symbols are used:



General warning symbol



Warning - dangerous high voltage



Warning - hot surfaces



Warning - automatic start-up of the unit



Biohazard warning

The warnings are structured as follows:

SIGNAL WORD

Description of the type and source of danger

Here you will find the possible consequences of ignoring the warning

> Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

DANGER

Immediate danger of severe injury or death

WARNING

Possible danger of severe injury or death

CAUTION

Risk of minor injuries

NOTICE

Risk of extensive material/property damage

Other symbols

These symbols are used in the document and on or in the unit:



Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.



Refer to the accompanying electronic documents.



Product for the disinfection and cleaning of suction units



Wear protective gloves.



Wear protective goggles.



Use a face mask.



Disconnect all power from the unit.



Lower and upper temperature limits



Lower and upper humidity limits



Protective ground connection



CE labelling with the number of the notified body



Conformity mark for the United Kingdom of Great Britain and Northern Ireland, with the number of the designated authority



Ukrainian conformity mark with registration number



Serial number



Order number



Medical device



Health Industry Bar Code (HIBC)



Manufacturer

1.2 Copyright information

All circuits, processes, names, software programs and units mentioned in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from the copyright owner.

Safety 2

The unit has been developed and designed in such a way that dangers are effectively ruled out if used in accordance with the Intended Use. Despite this, the following residual risks can remain:

- Personal injury due to incorrect use/misuse
- Personal injury due to mechanical effects
- Personal injury due to electrical shock
- Personal injury due to radiation
- Personal injury due to fire
- Personal injury due to thermal effects on skin
- Personal injury due to lack of hygiene, e.g. infection

2.1 Intended purpose

The suction unit provides the dental treatment unit with vacuum and volume flow.

22 Intended use

Working in combination with the suction unit with treatment unit, suction handpiece and cannula, the media used in dental treatment (e.g. water. saliva, dentine and amalgam) are removed by suction for disposal.

This unit is technically suitable for the aspiration of nitrous oxide (laughing gas). However, when assembling a system for aspiration of nitrous oxide, it is important to ensure that the other components in the system are also suitable for this purpose. Those responsible for setting up the system must assess this and approve and release the system for the aspiration of nitrous oxide.



Operation with nitrous oxide is only permitted if the exhaust air is transported from the unit to the outside of the building.

2.3 Improper use

- Do not use this device to aspirate flammable or explosive mixtures.
- Do not use the unit as a vacuum cleaner.
- Do not use chemicals containing chlorine or foaming chemicals.
- Operation in operating theatres of explosive areas is not permissible.

Systems, connection with 2.4 other devices

Additional devices connected with medical electrical devices must be proven to conform with their corresponding IEC or ISO standards. All configurations must continue to comply with the standard requirements for medical systems (see IEC 60601-1).

Whoever connects additional devices to medical electrical devices automatically becomes the system configurator and is responsible for ensuring that the system corresponds with the standard requirements for systems. Local laws take priority over the requirements outlined above.

2.5 General safety information

- Always comply with the specifications of all guidelines, laws, and other rules and regulations applicable at the site of operation for the operation of this unit.
- Check the function and condition of the unit prior to every use.
- Do not convert or modify the unit.
- Comply with the specifications of the Installation and Operating Instructions.
- The Installation and Operating Instructions must be accessible to all operators of the unit at all times.

2.6 Specialist personnel

Operation

Unit operating personnel must ensure safe and correct handling based on their training and knowledge.

 Instruct or have every operator instructed in handling the unit.

Installation and repairs

- The manufacturer recommends that installation, readjustments, alterations, upgrades and repairs be carried out either by the manufacturer itself or by a qualified specialist authorised by the manufacturer.

2.7 Notification requirement of serious incidents

The operator/patient is required to report any serious incident that occurs in connection with the device to the manufacturer and to the competent authority of the Member State in which the operator and/or patient is established/resident.

28 Electrical safety

- Comply with all the relevant electrical safety regulations when working on the unit.
- Never touch the patient and unshielded plug connections on the unit at the same time.
- Replace any damaged cables or plugs immediately.

Observe the EMC rules concerning medical devices

- The unit is intended for use in professional healthcare facilities (in accordance with IEC 60601-1-2). If the appliance is operated in another environment, potential effects on electromagnetic compatibility must be taken into account.
- Do not operate the unit in the vicinity of HF surgical instruments or MRT equipment.
- Maintain a minimum distance of at least 30 cm between the unit and other electronic devices.
- Note that cable lengths and cable extensions have effects on electromagnetic compatibility.

No maintenance measures are required to maintain the EMC basic safety.



NOTICE

Negative effects on the EMC due to non-authorised accessories

- > Use only those accessories named or approved by the manufacturer.
- Using any other accessories may result in increased electromagnetic interference emissions or the unit having reduced electromagnetic immunity, leading to an erroneous operation mode.





NOTICE

Erroneous operation mode due to use immediately adjacent to other devices or with other stacked devices

- Do not stack the unit together with other devices.
- If this is unavoidable, the unit and other devices should be monitored in order to ensure that they are working correctly.



NOTICE

Reduced performance characteristics due to insufficient distance between unit and portable HF communication devices

Keep a distance of at least 30 cm between the unit (including parts and cables of the unit) and portable HF communication devices (wireless units) (including their accessories such as antenna cables and external antennas).

- Decontaminate potentially contaminated parts before disposing of them.
- Uncontaminated parts (e.g. electronics, plastic and metal parts etc.) should be disposed of in accordance with the local waste disposal regulations.
- If you have any questions about the correct disposal of parts, please contact your dental trade supplier.



An overview of the waste keys for Dürr Dental products can be found in the download area:



http://gr.duerrdental.com/P007100155

2.9 Only use original parts

- Only use accessories and optional articles named or authorised by the manufacturer.
- Only use only original wear parts and replacement parts.

2.10 Transport

The original packaging provides optimum protection for the unit during transportation. If required, the original packaging for the unit can be ordered.



The manufacturer and the distributor do not accept liability, even during the warranty period, for damage during transportation due to improper packaging.

- Only transport the unit in its original packaging.
- Keep the packing materials out of the reach of children.

2.11 Disposal



The unit may be contaminated. Instruct the company disposing of the waste to take the relevant safety precautions.

Product description

Overview




Fig. 1: V 300 S

- 1 Suction unit
- 2 Condensation separator
- 3 Suction hose
- 4 Exhaust air hose (aluminium)
- 5 Waste water hose LW 20
- 6 Set of connection fittings
- 7 Mains cable




Fig. 2: VS 300 S

- 1 Combination suction unit
- 2 Suction hose
- 3 Exhaust air hose (aluminium)
- 4 Waste water hose LW 20
- 5 Set of connection fittings
- 6 Mains cable

3.1 Scope of delivery

The following items are included in the scope of delivery (possible variant-specific deviations due to country-specific requirements and/or import regulations):

V 300 S

V 300 S, 230 V,	1~, 50 Hz	7119-01
V 300 S, 230 V,	1~, 50/60 Hz	7119-02

- Connector set
- Mains cable

V 300 S. 230 V. 1~, 50 Hz..... 7119-01/002 V 300 S, 230 V, 1~, 50/60 Hz.... 7119-02/002

Set of connection fittings

- Suction hose LW 30, grey
- Exhaust air hose LW 30, aluminium
- Waste water hose LW 20
- Condensation separator
- Connector set
- Mains cable

VS 300 S

VS 300 S, 230 V,	1~, 50 Hz	7122-01
VS 300 S, 230 V,	1~, 50/60 Hz	7122-02

- Connector set
- Mains cable

VS 300 S, 230 V, 1~, 50 Hz, with rinsing unit 7122-04

- Protective strainer with connecting parts
- Connector set
- Mains cable

VS 300 S, 230 V, 1~, 50 Hz 7122-01/002 VS 300 S, 230 V, 1~, 50/60 Hz.. 7122-02/002 VS 300 S, 230 V, 1~, 50 Hz, with rinsing unit 7122-04/002

- Set of connection fittings
- Connector set
- Suction hose LW 30, grey
- Exhaust air hose LW 30, aluminium
- Waste water hose LW 20
- OroCup
- Mains cable

3.2 Optional items

The following optional items can be used with the device:

V 300 S

8 I EN

Wall bracket 7130-190-00

Noise reduction hood	7122200000
Condensation separator kit for housing	7119-701-20
Bacteria filter with accessories	
Bacteria filter with housing	
Ventilation kit for cabinet installa-	
tion	7122-981-51
Console for floor-mounted installa-	
tion	7130-191-00
VS 300 S	
Wall bracket	7130-190-00
Noise reduction hood	7100000000
. 10.00 . 0000	1122200000
Bacteria filter with accessories	
	7120-143-00
Bacteria filter with accessories	7120-143-00
Bacteria filter with accessories Bacteria filter with housing Rinsing unit conversion set for VS 300 S and VSA 300 S	7120-143-00 7120100000 7100-120-53
Bacteria filter with accessories Bacteria filter with housing Rinsing unit conversion set for	7120-143-00 7120100000 7100-120-53
Bacteria filter with accessories Bacteria filter with housing Rinsing unit conversion set for VS 300 S and VSA 300 S Rinsing unit II	7120-143-00 7120100000 7100-120-53 7100-250-50
Bacteria filter with accessories Bacteria filter with housing Rinsing unit conversion set for VS 300 S and VSA 300 S Rinsing unit II	7120-143-00 7120100000 7100-120-53 7100-250-50
Bacteria filter with accessories Bacteria filter with housing Rinsing unit conversion set for VS 300 S and VSA 300 S Rinsing unit II	7120-143-00 7120100000 7100-120-53 7100-250-50 7122-981-51

33 Consumables

The following materials are consumed during operation of the device and must be ordered separately:

Orotol plus (2.5 litre bottle) CDS110P6150 Orotol plus pH 7 (2.5-litre bottle). CDS117A6150 MD 555 cleaner (2.5 litre bottle). CCS555C6150 MD 555 cleaner organic (2.5-litre

3.4 Wear parts and replacement parts

The following working parts need to be changed at regular intervals (refer to the "Maintenance" section):

Nonreturn valve (pack of 3) 7128-100-03E



Information about replacement parts is available from the portal for authorised specialist dealers at: www.duerrdental.net

Technical data 4

4.1 V 300 S

Electrical data		7119-01	7119	9-02	
Rated voltage	V	230, 1~	230	, 1~	
Mains frequency	Hz	50	50	60	
Nominal current	А	2.9	2.9	3.7	
Starting current, approx.	А	~ 12	~ 12	~ 12	
Motor protection		Motor winding	overheat prot (±5 °C)	ector 160 °C	
Rated power	W	580	580	800	
Type of protection			IP 20		
Protection class			I		
Protective low voltage	V		24 ~		
Output	VA		4		
Connections					
Suction connection, DürrConnect special	mm		Ø 30		
Exhaust air connection (external)	mm		Ø 30		

Media				
Max. number of users			1	
Max. flow rate with unimpeded flow	l/min	740	740	820
Max. suction system pressure *	mbar/hPa		-200	

Depending on unit type

General data				
Max. number of users			1	
Duty cycle	%		100	
Heat generation rate	kJ/h	2088	2088	2880
Dimensions (H x W x D) *	cm		38 x 31 x 32	
Weight, approx. without housing with housing	kg kg		12.4 20.4	
Noise level ** approx. without housing with housing	dB(A) dB(A)	60 50	60 50	62 53

Values without accessories and add-on parts

Sound pressure level in accordance with ISO 3746

Ambient conditions during storage an	d transport	
Temperature	°C	-10 to +60

Ambient conditions during storage an	d transport	
Relative humidity	%	< 95
Ambient conditions during operation		
Temperature	°C	+10 to +40
Relative humidity	%	< 70
Classification		
Medical Device Class (MDR)		lla
Electromagnetic compatibility (EMC) Interference emission measurements		
High-frequency emissions in accordance	with CISPR 11	Group 1 Class B
Interference voltage at the power supply CISPR 11:2009+A1:2010	connection	Compliant
Electromagnetic interference radiation CISPR 11:2009+A1:2010		Compliant
Emission of harmonics IEC 61000-3-2:2005+A1:2008+A2:2009)	Compliant
Voltage changes, voltage fluctuations and sions IEC 61000-3-3:2013	d flicker emis-	Compliant
Electromagnetic compatibility (EMC) Interference immunity measurements		
Immunity to electrostatic discharge IEC 61000-4-2:2008		Compliant
Immunity to high-frequency electromagn IEC 61000-4-3:2006+A1:2007+A2:2010		Compliant
Immunity to near fields of wireless HF co devices IEC 61000-4-3:2006+A1:2007+A2:2010		Compliant
Immunity to fast electrical transients/burs voltage IEC 61000-4-4:2012	sts – AC mains	Compliant
Immunity to electrical fast transients/burs SIP/SOP ports IEC 61000-4-4:2012	sts – I/O,	Compliant
Immunity to interference, surges IEC 61000-4-5:2005		Compliant
Immunity to conducted disturbances, inc frequency fields – AC mains voltage IEC 61000-4-6:2013	duced by radio-	Compliant
Immunity to conducted disturbances, inc frequency fields – SIP/SOP ports IEC 61000-4-6:2013	duced by radio-	Compliant



Electromagnetic compatibility (EMC) Interference immunity measurements Immunity to power frequency magnetic fields Compliant IEC 61000-4-8:2009 Immunity to voltage dips, short interruptions and voltage Compliant variations IEC 61000-4-11:2004

Radio service Frequency band MHz Test level V/m TETRA 400 380 - 390 27 GMRS 460 FRS 460 FRS 460 FRS 460 FRS 460 FRS 460 FRS 460 430 - 470 28 LTE band 13, 17 704 - 787 9 GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5 800 - 960 28 CDMA 1900 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS 1700 - 1990 28 Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7 2400 - 2570 28 WLAN 802.11 a/n 5100 - 5800 9	Immunity to interference table, near fields of wireles	s HF communication devi	ces
GMRS 460 FRS 460 430 - 470 28 LTE band 13, 17 704 - 787 9 GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5 GSM 1800 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	Radio service		
FRS 460 430 - 470 28 LTE band 13, 17 704 - 787 9 GSM 800/900 TETRA 800 iDEN 820 800 - 960 28 CDMA 850 LTE band 5 GSM 1800 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	TETRA 400	380 - 390	27
GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5 GSM 1800 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7		430 - 470	28
TETRA 800 iDEN 820 CDMA 850 LTE band 5 GSM 1800 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	LTE band 13, 17	704 - 787	9
CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	TETRA 800 iDEN 820 CDMA 850	800 - 960	28
WLAN 802.11 b/g/n 2400 - 2570 28 RFID 2450 LTE band 7	CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25	1700 - 1990	28
WLAN 802.11 a/n 5100 - 5800 9	WLAN 802.11 b/g/n RFID 2450	2400 - 2570	28
	WLAN 802.11 a/n	5100 - 5800	9

Electromagnetic compatibility (EMC) Interference immunity measurements on the supply input	
Immunity to fast electrical transients/bursts – AC mains voltage IEC 61000-4-4:2012 ± 2 kV 100 kHz repetition rate	Compliant
Immunity to surges, line-to-line IEC 61000-4-5:2005 ± 0.5 kV, ± 1 kV	Compliant
Immunity to surges, line-earth IEC 61000-4-5:2005 ± 0.5 kV, ± 1 kV, ± 2 kV	Compliant

Electromagnetic compatibility (EMC)

Interference immunity measurements on the supply input

Immunity to conducted disturbances, induced by radio-

frequency fields - AC mains voltage

IEC 61000-4-6:2013

3 V

0.15-80 MHz Compliant

6 V

ISM frequency bands 0.15-80 MHz

80% AM at 1 kHz

Immunity to voltage dips, short interruptions and voltage

variations

IEC 61000-4-11:2004

Compliant

Electromagnetic compatibility (EMC)

Interference immunity measurements SIP/SOP

Immunity to electrostatic discharge

IEC 61000-4-2:2008

± 8 kV contact

 \pm 2kV, \pm 4 kV, \pm 8 kV, \pm 15 kV air

Immunity to electrical fast transients/bursts - I/O,

SIP/SOP ports

IEC 61000-4-4:2012 + 1 kV

100 kHz repetition rate

Immunity to impulse voltages, conductor to earth

IEC 61000-4-5:2005

 $\pm 2 kV$

Immunity to conducted disturbances, induced by radio-

frequency fields - SIP/SOP ports

IEC 61000-4-6:2013

3 V 0.15-80 MHz

6 V

ISM frequency bands 0.15-80 MHz

80% AM at 1 kHz

Compliant

Compliant

Compliant

4

50

4.2 **VS 300 S**

Electrical data		7122-01 7122-04	7122	2-02
Rated voltage	V	230, 1~	230	, 1~
Mains frequency	Hz	50	50	60
Nominal current	Α	2.9	2.9	3.7
Starting current, approx.	А	~ 12	~ 12	~ 12
Motor protection		Motor winding	overheat prot (±5 °C)	ector 160 °C
Rated power	W	580	580	800
Type of protection			IP 20	
Protection class			I	
Protective low voltage	V		24 ~	
Output	VA		4	
Connections				
Suction connection, DürrConnect special	mm		Ø 30	
Exhaust air connection (external)	mm		Ø 30	
Drain connection, DürrConnect	mm		Ø 20	
Media				
Max. flow rate with unimpeded flow	l/min	710	710	800
Max. suction system pressure *	mbar/hPa		-200	

*	Depending on unit type

Max. rate of flow of fluids

Max. suction height

General data				
Max. number of users			1	
Duty cycle	%		100	
Heat generation rate	kJ/h	2088	2088	2880
Dimensions (H x W x D) *	cm		38 x 31 x 32	
Weight, approx. without housing with housing	kg kg		12.6 20.6	
Noise level ** approx. without housing with housing	dB(A) dB(A)	60 50	60 50	62 53

I/min

cm

Values without accessories and add-on parts

Sound pressure level in accordance with ISO 3746

Ambient conditions during storage and	transport	
Temperature	°C	-10 to +60
Relative humidity	%	< 95
Ambient conditions during operation		
Temperature	°C	+10 to +40
Relative humidity	%	< 70
Classification		
Medical Device Class (MDR)		lla
Electromagnetic compatibility (EMC) Interference emission measurements		
High-frequency emissions in accordance v	with CISPR 11	Group 1 Class B
Interference voltage at the power supply c CISPR 11:2009+A1:2010	connection	Compliant
Electromagnetic interference radiation CISPR 11:2009+A1:2010		Compliant
Emission of harmonics IEC 61000-3-2:2005+A1:2008+A2:2009		Compliant
Voltage changes, voltage fluctuations and flicker emissions IEC 61000-3-3:2013		Compliant
Electromagnetic compatibility (EMC) Interference immunity measurements		
Immunity to electrostatic discharge IEC 61000-4-2:2008		Compliant
Immunity to high-frequency electromagnet IEC 61000-4-3:2006+A1:2007+A2:2010	tic fields	Compliant
Immunity to near fields of wireless HF comdevices IEC 61000-4-3:2006+A1:2007+A2:2010	nmunication	Compliant
Immunity to fast electrical transients/bursts voltage IEC 61000-4-4:2012	s – AC mains	Compliant
Immunity to electrical fast transients/bursts SIP/SOP ports IEC 61000-4-4:2012	s – I/O,	Compliant
Immunity to interference, surges IEC 61000-4-5:2005		Compliant
Immunity to conducted disturbances, indufrequency fields – AC mains voltage IEC 61000-4-6:2013	uced by radio-	Compliant



Electromagnetic compatibility (EMC) Interference immunity measurements	
Immunity to conducted disturbances, induced by radio- frequency fields – SIP/SOP ports IEC 61000-4-6:2013	Compliant
Immunity to power frequency magnetic fields IEC 61000-4-8:2009	Compliant
Immunity to voltage dips, short interruptions and voltage variations IEC 61000-4-11:2004	Compliant

Immunity to interference table, near fields of wireless HF communication devices		
Radio service	Frequency band MHz	Test level V/m
TETRA 400	380 - 390	27
GMRS 460 FRS 460	430 - 470	28
LTE band 13, 17	704 - 787	9
GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE band 5	800 - 960	28
GSM 1800 CDMA 1900 GSM 1900 DECT LTE band 1, 3, 4, 25 UMTS	1700 - 1990	28
Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE band 7	2400 - 2570	28
WLAN 802.11 a/n	5100 - 5800	9

Electromagnetic compatibility (EMC) Interference immunity measurements on the supply input	
Immunity to fast electrical transients/bursts – AC mains voltage IEC 61000-4-4:2012 ± 2 kV 100 kHz repetition rate	Compliant
Immunity to surges, line-to-line IEC 61000-4-5:2005 ± 0.5 kV, ± 1 kV	Compliant
Immunity to surges, line-earth IEC 61000-4-5:2005 ± 0.5 kV, ± 1 kV, ± 2 kV	Compliant

Electromagnetic compatibility (EMC)

Interference immunity measurements on the supply input

Immunity to conducted disturbances, induced by radio-

frequency fields - AC mains voltage

IEC 61000-4-6:2013

3 V

0.15-80 MHz

Compliant

6 V

ISM frequency bands

0.15-80 MHz 80% AM at 1 kHz

Immunity to voltage dips, short interruptions and voltage

variations

IEC 61000-4-11:2004

Compliant

Electromagnetic compatibility (EMC)

Interference immunity measurements SIP/SOP

Immunity to electrostatic discharge

IEC 61000-4-2:2008

± 8 kV contact

 \pm 2kV, \pm 4 kV, \pm 8 kV, \pm 15 kV air

Immunity to electrical fast transients/bursts - I/O,

SIP/SOP ports

IEC 61000-4-4:2012

+ 1 kV

100 kHz repetition rate

Immunity to impulse voltages, conductor to earth

IEC 61000-4-5:2005

 $\pm 2 kV$

Immunity to conducted disturbances, induced by radio-

frequency fields - SIP/SOP ports

IEC 61000-4-6:2013

3 V

0.15-80 MHz

6 V

ISM frequency bands

0.15-80 MHz

80% AM at 1 kHz


Compliant

Compliant

Compliant

4.3 Type plate

The type plate is is located on the noise reduction hood.



1 Type plate


Evaluation of conformity 4.4

This device has been subjected to conformity acceptance testing in accordance with the current relevant European Union guidelines. This equipment conforms to all relevant requirements.



5 Function

5.1 V 300 S



- 1 Motor
- 2 Auxiliary air nozzle
- 3 Exhaust air muffler
- 4 Mains connection (power plug)
- 5 Control electronics
- 6 Control connection
- 7 Control signal output
- 8 Suction connection
- 9 Turbine wheel
- 10 Exhaust air connection

When treating their patients with dental instruments, dentists, dental assistants and hygienists all face increased exposure to infectious aerosols and droplets. Studies have shown that intraoral spray mist extraction can effectively eliminate these aerosols and droplets. According to one study, a flow rate of at least 300 l/min reduces particle emissions when using a turbine to below the limit of detection (study by Koch & Graetz).




The V-suction unit is suitable for use in dry air suction systems. The advantage of this system is that the suction units can be installed in any suitable room, regardless of the actual connection layout and routeing of the lines. The necessary air flow and vacuum are generated by a rapidly rotating impeller. An auxiliary air nozzle on the turbine housing protects the suction units against overheating.

When an appropriate vacuum for the machine is applied, approx. 300 l/min of air is sucked in through the suction cannula.

On the vacuum side the V-suction unit is equipped with a condensation separator that collects any condensation arising within the pipe system and transports it away to the outside.

The exhaust air from the suction unit should be guided out of the building (via the roof where possible). We recommended the installation of a bacteria filter in the exhaust air line. In addition, it is possible to install a noise-reducing muffler in the exhaust air line in order to reduce the amount of noise generated by the unit and by the air flow.

5.2 VS 300 S



- 1 Motor
- 2 Auxiliary air nozzle
- 3 Exhaust air muffler
- 4 Mains connection (power plug)
- 5 Control electronics
- 6 Control connection
- 7 Control signal output
- 8 Waste water connection
- 9 Separation
- 10 Suction connection
- 11 Turbine wheel
- 12 Exhaust air connection

When treating their patients with dental instruments, dentists, dental assistants and hygienists all face increased exposure to infectious aerosols and droplets. Studies have shown that intraoral spray mist extraction can effectively eliminate these aerosols and droplets. According to one study, a flow rate of at least 300 l/min reduces particle emissions when using a turbine to below the limit of detection (study by Koch & Graetz).



The VS suction units are used in wet suction systems. The suction units can be installed on the same floor as the treatment units or on the floor underneath. The necessary air flow and vacuum are generated by a rapidly rotating impeller. An auxiliary air nozzle on the turbine housing protects the suction units against overheating.

When an appropriate vacuum for the machine is applied, approx. 300 l/min of air is sucked in through the suction cannula.

The impeller, the separation turbine and the waste water pump are driven by the motor.

The mixture of liquids, solid particles and air drawn in passes through the inlet connection and into the suction unit. The coarse filter holds back the solid particles.

Inside the separation unit, the aspirated fluids and solid particles pass through a two-stage separation system and are separated from the suction air. This separation system consists of a cyclone separator and a separation turbine. The suction process runs continuously.

The aspirated mixture flows into the cyclonic separator, where it is set into a spiral motion. In this first stage, the resulting centrifugal forces force the fluid constituents and any remaining solid particles against the outside wall of the separation chamber of the cyclone separator. This initially only effects a "coarse separation" of the fluid. In the subsequent second stage, the separation turbine effects "fine separation" of the remaining liquid from the air flow that has carried it to this point.

The waste water pump transports the liquid from the centrifuge together with the fine solid particles through the waste water connection into the central waste water network. A diaphragm valve is located in the waste water connection to prevent fluid from the drain being sucked back in.

The exhaust air from the suction unit should be guided out of the building (via the roof where possible). We recommended the installation of a bacteria filter in the exhaust air line. In addition, it is possible to install a noise-reducing muffler in the exhaust air line in order to reduce the amount of noise generated by the unit and by the air flow.



6 Requirements

Depending on the suction system, different installation options are available.



Further information can be found in our suction planning information leaflet.

Order number9000-617-03/..

6.1 Installation/setup room

The room chosen for set up must fulfil the following requirements:

- Closed, dry, well-ventilated room
- Should not be a room made for another purpose (e. g. boiler room or wet cell)
- When installing in a cabinet the inlet and outlet ventilation slots must be present; minimum free cross-section at least 120 cm².
- Forced ventilation (fan) must be provided if there is a risk that the recommended room air temperature could be exceeded. The air flow performance must be at least 2 m³/min.
- Do not cover cooling slots or openings with housing installations; ensure sufficient clearance to the openings to permit sufficient cooling.
- Mains cable plug connections must be freely accessible so they can be quickly disconnected if there is any danger.

6.2 Setup options

The following options for setting up the unit are available:

- Wall installation using a Dürr Dental wall mounting
- In a ventilated cabinet
- In a Dürr Dental noise reducing housing

6.3 Pipe materials

Only use vacuum-sealed HT-waste pipes manufactured from the following materials:

- Polypropylene (PP),
- Chlorinated polyvinyl chloride (PVC-C),
- Plasticizer-free polyvinyl chloride (PVC-U),
- Polyethylene (PE).

The following materials must not be used:

- Acrylonitrile-butadiene-styrene (ABS),
- Styrene copolymer blends (e.g. SAN + PVC).

6.4 Hose materials

For waste connections and suction lines only use the following hose types:

- Flexible spiral hoses made of PVC with integrated spiral or equivalent hoses
- Hoses that are resistant to dental disinfectants and chemicals



Plastic hoses will display signs of ageing over time. Therefore, they should be inspected regularly and replaced as necessary.

The following types of hoses must not be used:

- Rubber hoses
- Hoses made completely of PVC
- Hoses that are not sufficiently flexible

6.5 Information about electrical connections

- Ensure that the electrical connections to the mains power supply are established in accordance with current valid national and local regulations and standards governing the installation of low voltage units in medical facilities.
- Observe the current consumption of the devices that are to be connected.

Electrical fusing

LS switch 16 A, characteristic B, C and D in accordance with 60898.

6.6 Information about connecting cables

Mains supply cable

Only use the supplied mains cable to connect the device.

Control cable

24 V protective low voltage for:

- Hose manifold
- Place selection valve
- Spittoon valve

Installation type	Line layout (minimum requirements)
Fixed installation	 Shielded sheathed cable (e.g. (N)YM (St)-J)

Installation type	Line layout (minimum requirements)
Flexible	 PVC data cable with shielded cable sheath- ing, as used for tele- communications and IT processing systems (e.g. type LiYCY)
	or - Lightweight PVC control cable with shielded cable sheathing



Connect the shielding of the cables in accordance with the regulations.

System components

The system components listed below are required or recommended for various procedures or for installation.

7.1 OroCup



The OroCup care system is a closed dosing system for simple preparation and suctioning of disinfectants and special cleaners.

Use OroCup to clean and disinfect the suction system with all of its components and the spittoon.

7.2 Condensate separator

The condensation separator allows targeted and reliable extraction of condensation in the suction line system.

The condensation separator is installed at the lowest point of the suction line system.

7.3 Rinsing unit

It is recommended that the suction system is equipped with a rinsing unit, e.g. in the treatment unit. The rinsing unit provides a small amount of water during aspiration. This dilutes the aspirated fluids (blood, saliva, rinsing water etc.), which can then be transported more effectively.

7.4 Flow accelerator

In order to keep the suction system free of deposits, a flow accelerator can be fitted in conjunction with a spittoon valve. When using a bowl rinse system, water will collect before the flow accelerator. The next time suction takes place using the large cannula, the collected fluid is transported in surges and at high speed to the suction system. This ensures automatic cleaning of the suction pipes.

7.5 Amalgam separator


The amalgam separator is designed to separate out and trap the heavy metal particles and amalgam dust that the suction unit aspirates from drilled fillings. The amalgam separator is installed in

the drain behind the separation unit of the suction unit. The amount of fluid coming from the suction unit must not exceed the maximum permitted quantity of fluid that can be handled by the amalgam separator. Depending on the installation and on national regulations, a second amalgam separator may need to be installed.

7.6 Surge tank

If the suction unit is combined with an amalgam separator, this requires the installation of a surge tank. The surge tank reduces pressure peaks caused by the waste water pump of the suction unit and acts as a buffer against temporary rises in the volume of water.

The surge tank can also be used if the waste water is fed directly into the building waste water system. this case the waste water from the suction unit is diverted to the building drainage system under zero pressure.



- 1 CA₄1
- 2 Pressure equalisation tank
- 3 VS 300 S combination suction unit

7.7 Bacteria filter

For hygienic reasons, we recommend the installation of a bacteria filter in the exhaust air line. If the unit is installed in the surgery and the exhaust air cannot be discharged to the outdoors, it is essential to install a bacteria filter. Depending on the type and condition of the bacteria filter, it will need to be replaced every 1-2 years at the latest.



The separation integrated in the system does not retain bacteria; this is why we recommend installing a suitable filter in the exhaust air system.

7.8 Noise reduction

If the noise level from the exhaust air vent or the flow noise generated is too high, noise reduction can be installed in the exhaust air line.


8 Installation



The actual connection can vary depending on the chosen installation option. The connection shown is only an example.

8.1 Mounting the condensation separator on the V 300 S


- 1. Use two bolts and two nuts to screw the condensation separator onto the retaining plate.
- 2. Use two bolts to screw the retaining plate with the condensation separator onto the noise reduction hood.
- 3. Push the hose sleeve onto the suction connection on the suction unit and secure it with the securing ring.
- 4. Push the connecting hose onto the hose adapters on the suction unit and on the condensation separator.



8.2 Installation and routeing of hoses and pipes


- Establish connections between the pipe system and the unit using the flexible hoses supplied. This will prevent vibrations from being transmitted to the pipe system.
- The connection between the pipe line and unit suction connection should be kept as short as possible and straight, without bends.
- Install the drain hoses with a downward gradient so that the waste water can drain off.
- Waste water connections must be implemented in accordance with applicable local and national regulations.

V 300 S




- 1 O-ring Ø 30x2 mmr
- 2 Hose connector Ø 30 mm
- 3 Hose clip 25-40 mm
- 4 Waste air pipe (aluminium)Ø 30 mm inside
- 5 Suction hose Ø 30 mm (internal)
- 6 Hose sleeve
- 7 Hose clip Ø 28 mm
- 8 Waste water hose Ø 20 mm (internal)
- 9 Hose sleeve Ø 20 mm
- 10 Securing ring
- 11 O-ring Ø 20x2 mm
- 12 Connector Ø 36 mm (external)

VS 300 S




- O-ring Ø 30x2 mmr 1
- 2 Hose connector Ø 30 mm
- 3 Hose clip 25-40 mm
- 4 Waste air pipe (aluminium)Ø 30 mm inside
- 90° elbow, DN 30 5
- 6 Securing ring
- 7 Suction hose Ø 30 mm (internal)
- 8 Connector Ø 36 mm (external)
- 9 O-ring Ø 20x2 mm
- 10 Securing ring
- 11 Hose sleeve Ø 20 mm
- 12 Hose sleeve
- 13 Hose clip Ø 28 mm
- 14 Waste water hose Ø 20 mm (internal)
- 15 Hose sleeve Ø 30 mm
- 16 45° elbow, DN 30

8.3 Rinsing unit water connections




- 1. Slide the clamping ring approx. 1.5 cm down the water hose.
- Dürr Dental recommends a water hose with an internal diameter of 2 mm; material: TPU, 87 Shore A, test certificate in accordance with German KTW ("contact with potable water") guideline.
- Plug the water hose onto the water connector.
- 3. Use a suitable tool to slide the clamping ring to just before the end of the water hose.



- Attach a T-piece for the water hose with a diameter of Ø 4 mm or Ø 6 mm in the water supply.
- 5. Attach the water hose with the sleeve piece, clamping ring and union nut to the T-piece.

 Alternatively, attach the water hose with the adapter part, seal, R3/4" screw connection, sleeve piece, double-tapered ring and union nut to a water tap.




Electrical connections 9

NOTICE

Short circuit due to defective lead

Do not route wires near hot surfaces.

- 1. Before connecting, check that the power supply voltage matches the voltage specifications on the type plate.
- 2. Connect control line to control connection.
- 3. Plug the mains cable into the device.



- Power plug
- X1 Mains connection on the circuit board
- X2 Motor connection
- X3 Control connection 24 V AC / max. 80 mA
- X4 Control signal output 24 V AC / max. 20 mA

10 Commissioning



- In many countries technical medical products and electrical devices are subject to regular checks at set intervals. The owner must be instructed accordingly.
- 1. Turn on the unit power switch or the main surgery switch.
- 2. Carry out a function check of the system.
- Check all connections for leak tightness.
- 4. Carry out an electrical safety check in accordance with applicable regulations (e.g. regulations concerning set up, operation and application of medical devices) and record the results as appropriate, e.g. in the technical log book.
- 5. Carry out and document the instruction and handover for the unit.



A sample handover report is included in the attachment.



Usage



When using prophy powders, water-soluble Lunos Prophy Powders are recommended in order to protect the suction systems (Dürr Dental).

Disinfection and cleaning

The following tasks are required for the disinfection and cleaning of the suction system:

"11.1 Suctioning water"	After every treatment
"11.2 Disinfection"	Daily in the evening after the end of treatment, With increased workloads, before the midday break and after the end of treatment
"11.3 Cleaning"	Allow to act for at least 2 x per week before the start of treatment or during the midday break Alternatively, 5 x per week, with short reaction times



NOTICE

Device malfunctions or damage due to use of incorrect media

Guarantee claims may become invalid as a result.

- Do not use any foaming agents such as household cleaning agents or instrument disinfec-
- > Do not use abrasive cleaners.
- > Do not use agents containing chlorine.
- > Do not use any solvents like acetone.

The following should always be used:



- For disinfection and cleaning:
 - Orotol plus
 - Orotol plus pH 7
- For cleaning:
 - MD 555 cleaner
 - MD 555 cleaner organic

Only these products have been tested by Dürr Dental.



Comply with the instructions for use of the cleaning agent and disinfectant.

11.1 Suctioning water

After every treatment:

1. Suction up cold water (min. 0.5 litres) with the large and small suction hoses. Do this even if only the small suction hose was actually used during treatment.



Suction through the large suction hose causes a large amount of air to be drawn up, thereby considerably increasing the cleaning effect.



11.2 Disinfection

Disinfect and clean the suction system every evening after the end of treatment.



With increased workloads, clean and disinfect twice per day, e.g. before the midday break and after the end of treatment.

The following is required for the disinfection and cleaning of the suction system:



Non-foaming disinfectant/cleaning agent that is compatible with the materials

- Orotol plus
- Orotol plus pH 7



care system

e.g. OroCup



Disinfection and cleaning are described below with OroCup and Orotol plus.



Wear protective gloves.



Wear protective goggles.

1. To pre-clean, suction up 2 litres of water with the care system.



Using the solution in the care system:

Pour 2 sealing caps (40 ml) of Orotol plus into the OroCup.



Fill the OroCup with 2 litres of cold water.





Close the lid of the care system.



Mix the solution.



Using the solution:

1. Attach suction hoses to the care system and suction up 1 litre of solution.



Pour the remaining solution into the spittoon.



- Allow to work for a minimum of one hour or leave overnight.
- 4. When placing the system back into operation, suction up 2 litres of water.



The instructions are also available as a video:



11.3 Cleaning

Clean the suction system regularly:

- Allow to act for at least 2 x per week in the morning before the start of treatment or during the midday break, reaction time 30-120 minutes
- Alternatively: 5 x per week before the start of treatment, reaction time 5-15 minutes

The following is required for cleaning the suction system:



Special non-foaming detergent for suction systems that is compatible with the

- MD 555 cleaner
- MD 555 cleaner organic



care system

e.g. OroCup



Cleaning is described below with OroCup and MD 555 cleaner.



Wear protective gloves.



Wear protective goggles.

To pre-clean, suction up 2 litres of water with the care system.



Use the solution in the care system. Pour 5 caps (100 ml) of MD 555 cleaner into the care system.



3. Fill the OroCup with 2 litres of cold water.





4. Close the lid of the care system.



Mix the solution.



6. Attach suction hoses to the care system and suction up one litre of solution.



7. Pour the remaining solution into the spittoon.



8. Allow the solution time to act. For 2 x per week: 30–120 minutes For 5 x per week: 5-15 minutes

9. Suction up 2 litres of water after the reaction time has elapsed.



The instructions are also available as a video:





12 Maintenance



All maintenance work must be performed by a qualified expert or by one of our Service Technicians.



WARNING

Infection due to contaminated unit

- > Clean and disinfect the suction before working on the unit.
- > Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).



Prior to working on the unit or in case of danger, disconnect it from the mains.

12.1 V 300 S

Maintenance interval	Maintenance work
Every 1-2 years	Replace the exhaust air filter (where fitted). *
Every 2 years	Check the waste valve on the condensation separator for correct operation and replace it if necessary. *

* Only to be done by service technicians.

12.2 VS 300 S

Maintenance interval	Maintenance work
Every 4 weeks	Check the protective strainer on the suction connection of the unit and clean or replace it as required.
Annually	Check the non-return valve for correct operation and replace it if necessary. *
Every 1-2 years	> Replace the exhaust air filter (where fitted).
* Only to be performed by service technicians.	

12.3 Wear parts and replacement parts



Information about replacement parts is available from the portal for authorised specialist dealers at:

www.duerrdental.net


12.4 Cleaning the protective strainer



WARNING

Infection due to contaminated unit

- Clean and disinfect the suction before working on the unit.
- Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).
- Pull off the suction hose from the protective strainer.
- Pull off any hoses connected to the connection piece on the protective strainer.
- Pull out the protective strainer from the connection piece on the separation housing.
- 4. Clean the protective strainer.
- Push the protective strainer into the connection piece on the separation housing.
- 6. Reconnect all hoses that have been pulled off



1 Protective strainer

12.5 Intensive cleaning

Intensive cleaning can be carried out during the self-clean process of suction systems that are heavily contaminated with biofilms or if the suction power is greatly reduced.

Carry out this intensive cleaning at least twice a week.

Exposure time: 1 – 2 hours

The following is required for intensive cleaning of the suction system:



Special non-foaming detergent for suction systems that is compatible with the materials

- MD 555 cleaner
- MD 555 cleaner organic



care system
- e.g. OroCup



Cleaning is described below with OroCup and MD 555 cleaner.



Wear protective gloves.



Wear protective goggles.

1. To pre-clean, suction up 2 litres of water with the care system.



Use the solution in the care system. Pour 10 caps (200 ml) of MD 555 cleaner into the care system.



3. Fill the OroCup with 2 litres of cold water.





4. Close the lid of the care system.



5. Mix the solution.



6. Attach suction hoses to the care system and suction up one litre of solution.



7. Pour the remaining solution into the spittoon.



8. Allow the solution to act for 1 - 2 hours.

9. Suction up 2 litres of water after the reaction time has elapsed.



Troubleshooting

Tips for operators and service technicians



Any repairs exceeding routine maintenance may only be carried out by qualified personnel or our service.



WARNING

Infection due to contaminated unit

- > Clean and disinfect the suction before working on the unit.
- > Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).



Prior to working on the unit or in case of danger, disconnect it from the mains.

Error	Possible cause	Remedy
Device does not start	No mains voltage	 Check the mains voltage. * Check the fuses and replace if necessary. *
	Undervoltage	Measure the supply voltage; call an electrician if necessary.
	No start signal	Check the control voltage at the signal input. *
	Capacitor defective	Measure capacitance and replace if necessary. *
	Turbine is blocked by solid particles or sticky soiling	Disassemble the unit and clean the turbine. *
The unit generates unusual noises	Solid particles in the turbine chamber	Disassemble the unit and clean the turbine and housing.
Water leaking from the exhaust air connection	Membrane valve blocked	Check the membrane valve at the waste water connection and if necessary clean or replace. *
	Foam in turbine due to use of incorrect cleaning and disinfectant agents	Use non-foaming cleaning and disinfectant agents.
	Build-up of condensate in the exhaust air line	Check the pipe system; avoid over-cooling. *
	Waste water line/siphon trap blocked	Clean the waste water line/ siphon trap. *

Troubleshooting

Error	Possible cause	Remedy
Suction performance too low	Coarse filter blocked	Clean the coarse filter at the intake connection.
	Leak in the suction line	Check and if necessary establish leak-tightness of suction system and connec- tions. *
	Mechanical sluggishness of tur- bine caused by soiling	Disassemble the unit and clean the turbine. *

Only to be done by service technicians.

14 Transporting the unit



WARNING

Infection due to contaminated unit

- > Disinfect the unit before transport.
- Close all media connections.



Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

- Prior to disassembly, clean and disinfect the suction unit and the unit via aspiration of a suitable disinfectant approved by the manufacturer.
- Disinfect a defective unit using a suitable surface disinfection agent.
- 3. Seal all connections with sealing caps.
- **4.** Pack the unit securely in preparation for transport.




Fig. 3: V 300 S

- 1 Dummy bushing
- 2 Securing ring




Fig. 4: VS 300 S

- 1 Auxiliary air connection sealing cap
- 2 Rinse connection sealing cap
- 3 Dummy bushing
- 4 Water outflow sealing cap
- 5 Securing ring



15 Handover record

This document confirms that a qualified handover of the medical device has taken place and that appropriate instructions have been provided for it. This must be carried out by a qualified adviser for the medical device, who will instruct you in the proper handling and operation of the medical device.

Product name	Order number (REF)	Serial number (SN)
☐ Visual inspection of the p	oackaging for any damage	
 Unpacking the medical of 	device and checking for dam	age
□ Confirmation of the com	pleteness of the delivery	
Instruction in the proper instructions	handling and operation of the	e medical device based on the operating
Notoo		
Notes:		
Name of person receiving	instruction: Sign	ature:
Name and address of the	qualified advisor for the me	adical devices
ivallie and address of the c	qualified adviser for the me	dicai device.
Date of handover:		ature of the qualified adviser for the ical device:

40 | EN 7119100007L02 2501V004

16 Country representatives

Country

GB



Address

UK Responsible Person:

Duerr Dental (Products) UK Ltd. 14 Linnell Way Telford Way Industrial Estate Kettering, Northants NN 16 8PS

IJΑ



UA.TR.099

Уповноважений представник в Україні:

Приватне підприємство "Галіт" вул. 15 квітня, 6Є, с. Байківці,

Тернопільський р-н, 47711, Україна тел.: 0800 502 998; +38 050 338 10 64 www.galit.te.ua; e-mail: office@galit.te.ua

Виробник: Дюрр Дентал ЕсЕ Хьопфігхаймер Штрассе 17, Д-74321 Бітігхайм-Біссінген, Німеччина

email: info@duerrdental.com

CN

备案人/生产企业: DÜRR DENTAL SE 德国迪珥齿科股份公司

住所/生产地址: Höpfigheimer Str. 17, 74321 Bietigheim-Bissingen, Ger-

many 联系方式:

电话: +497142705-0 邮箱: info@duerrdental.com 网址: www.duerrdental.com

代理人/售后服务单位:迪珥医疗器械(上海)有限公司 住所:上海市长宁区天山路 641 号 2 号楼 (20 幢) 303 室

联系方式:

电话: +86 21 6381 0270 传真: +86 21 6381 0290 邮箱:info@duerr.cn

网址: http://www.duerrdental.com



Hersteller / Manufacturer:

DÜRR DENTAL SE Höpfigheimer Str. 17 74321 Bietigheim-Bissingen Germany

Fon: +49 7142 705-0 www.duerrdental.com info@duerrdental.com

