

Tyscor VS 2 Plus / VS 4



EN Installation and operating instructions

The current version of the installation and operating instructions is available in the Download Center:



<http://qr.duerdental.com/7188100021>

© 2019-2025 DÜRR DENTAL SE

Contents



Important information

1	About this document	3
1.1	Warnings and symbols	3
1.2	Copyright information	4
2	Safety	4
2.1	Intended purpose	4
2.2	Intended use	4
2.3	Improper use	4
2.4	General safety information	5
2.5	Combining devices safely	5
2.6	Specialist personnel	5
2.7	Notification requirement of serious incidents	5
2.8	Electrical safety	5
2.9	Only use original parts	5
2.10	Transport	6
2.11	Disposal	6
2.12	Cybersecurity measures (security context)	6



Product description

3	Overview	7
3.1	Scope of delivery	8
3.2	Optional items	8
3.3	Consumables	8
3.4	Wear parts and replacement parts	8
4	Technical data	9
4.1	Characteristic curves	11
4.2	Type plate	12
4.3	Evaluation of conformity	12
5	Function	13
5.1	Separation system	14
5.2	Radial blower	14
5.3	Bacteria filter (optional)	14
5.4	LEDs	14
5.5	Touch screen	14
5.6	Operating modes	14

5.7	Quick Start	14
5.8	Lag time	14
5.9	Eco Stop	14
5.10	VistaSoft Monitor (optional)	15
5.11	Grouping of multiple devices	15
5.12	OroCup	16



Assembly

6	Requirements	17
6.1	Installation/setup room	17
6.2	Setup options	17
6.3	Pipe materials	17
6.4	Hose materials	17
6.5	Information about electrical connections	17
6.6	Information about connecting cables	17
7	System components	18
7.1	Rinsing unit	18
7.2	Exhaust air filter / bacteria filter	18
7.3	Noise reduction	18
7.4	Surge tank	18
7.5	Flow accelerator	18
8	Installation	20
8.1	Installation and routeing of hoses and pipes	20
8.2	Installing the group of devices	21
8.3	Fitting the rinsing unit	21
8.4	Network connection	22
8.5	CAN bus	22
8.6	Electrical connections	22
8.7	PCB (main board) electrical connections	24
8.8	Testing suction power after installation	26
9	Commissioning	26
9.1	Setup wizard	26
9.2	Settings on the touch screen	26
9.3	Monitoring the unit via the network	29



Usage

10 LEDs	30
11 Operating the touch screen	31
11.1 Querying messages	31
12 Disinfection and cleaning	32
12.1 Suctioning water	33
12.2 Disinfection	33
12.3 Cleaning	34
13 Maintenance	36
13.1 Wear parts and replacement parts	36
13.2 Clean the filter at the device suction connection.	37
13.3 Confirm maintenance	37
13.4 Intensive cleaning	37



Troubleshooting

14 Tips for operators and service technicians	39
14.1 General faults	39
14.2 Fault with Tyscor Tandem Flex	42
14.3 Error messages	43
15 Transporting the unit	46



Appendix

16 Menu structure	47
17 Handover record	49
18 Country representatives	50

! Important information

1 About this document

These installation and operating instructions represent part of the unit.



The manufacturer and the distributor will not offer any guarantee or accept any liability for the safe operation and the safe functioning of the unit if the instructions and information in these installation and operating instructions are not complied with.

The German version of the installation and operating instructions is the original manual. All other languages are translations of the original manual. These installation and operating instructions apply to:

Tyscor VS 2 Plus

REF: 7188200100

Tyscor VS 4

REF: 7188100100

1.1 Warnings and symbols

Warnings

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The following warning symbols are used:



General warning symbol



Warning – dangerous high voltage



Warning - automatic start-up of the unit



Biohazard warning

The warnings are structured as follows:



SIGNAL WORD

Description of the type and source of danger

Here you will find the possible consequences of ignoring the warning

- Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

- **DANGER**
Immediate danger of severe injury or death
- **WARNING**
Possible danger of severe injury or death
- **CAUTION**
Risk of minor injuries
- **NOTICE**
Risk of extensive material/property damage

Other symbols

These symbols are used in the document and on or in the unit:



Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.



Product for the disinfection and cleaning of suction units



Wear protective gloves.



Disconnect all power from the unit.



Refer to the accompanying electronic documents.



Comply with the lower and upper temperature limits



Comply with the lower and upper humidity limits



CE labelling with the number of the notified body



Conformity mark for the United Kingdom of Great Britain and Northern Ireland, with the number of the designated authority



Authorised representative for Switzerland



Order number



Serial number



Medical device



Health Industry Bar Code (HIBC)



Manufacturer



Fragile, handle with care



This way up / store and transport in an upright position



Keep dry



Keep away from sunlight



Stacking limits

1.2 Copyright information

All circuits, processes, names, software programs and units mentioned in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from the copyright owner.

2 Safety

The unit has been developed and designed in such a way that dangers are effectively ruled out if used in accordance with the Intended Use. Despite this, the following residual risks can remain:

- Personal injury due to incorrect use/misuse
- Personal injury due to mechanical effects
- Personal injury due to electrical shock
- Personal injury due to radiation
- Personal injury due to fire
- Personal injury due to thermal effects on skin
- Personal injury due to lack of hygiene, e.g. infection

2.1 Intended purpose

The suction unit provides the dental treatment unit with vacuum and volume flow.

2.2 Intended use

Working in combination with the suction unit with treatment unit, suction handpiece and cannula, the media used in dental treatment (e.g. water, saliva, dentine and amalgam) are removed by suction for disposal.

This unit is technically suitable for the aspiration of nitrous oxide (laughing gas). However, when assembling a system for aspiration of nitrous oxide, it is important to ensure that the other components in the system are also suitable for this purpose. Those responsible for setting up the system must assess this and approve and release the system for the aspiration of nitrous oxide.



Operation with nitrous oxide is only permitted if the exhaust air is transported from the unit to the outside of the building.

2.3 Improper use

- Do not use this device to aspirate flammable or explosive mixtures.
- Do not use the unit as a vacuum cleaner.
- Do not use chemicals containing chlorine or foaming chemicals.
- Operation in operating theatres of explosive areas is not permissible.
- The suction unit must not be set up in the immediate surroundings of the patient (minimum distance: 1.5 m).

2.4 General safety information

- Always comply with the specifications of all guidelines, laws, and other rules and regulations applicable at the site of operation for the operation of this unit.
- Check the function and condition of the unit prior to every use.
- Do not convert or modify the unit.
- Comply with the specifications of the Installation and Operating Instructions.
- The Installation and Operating Instructions must be accessible to all operators of the unit at all times.

2.5 Combining devices safely

Take care when connecting units together or to parts of other systems as there is always an element of risk (e.g. due to leakage currents).

- Only connect units when there can be no question of danger to operator or to patient.
- Only connect units when it is safe to do so and when there is no risk of damage or harm to the surroundings.
- If it is not completely clear from the data sheet of the unit that such connections can be safely made or if you are in any doubt, always get a suitably qualified person (e.g. the relevant manufacturer) to verify that the setup is safe.

Where applicable, the requirements for medical products have been taken into account in the development and construction of the device. As a result, this device is suitable for installation within medical supply equipment.

- Where this device is integrated in other medical supply equipment, the requirements of European Union Medical Device Regulation 2017/745 and the relevant standards must be observed.

2.6 Specialist personnel

Operation

Unit operating personnel must ensure safe and correct handling based on their training and knowledge.

- Instruct or have every operator instructed in handling the unit.

The following groups are not permitted to operate or use a commercially operated unit:

- People without the necessary experience and knowledge
- People with reduced physical, sensory or mental capabilities
- Children

Installation and repairs

- The manufacturer recommends that installation, readjustments, alterations, upgrades and repairs be carried out either by the manufacturer itself or by a qualified specialist authorised by the manufacturer.

2.7 Notification requirement of serious incidents

The operator/patient is required to report any serious incident that occurs in connection with the device to the manufacturer and to the competent authority of the Member State in which the operator and/or patient is established/resident.

2.8 Electrical safety

- Comply with all the relevant electrical safety regulations when working on the unit.
- Never touch the patient and unshielded plug connections on the unit at the same time.
- Replace any damaged cables or plugs immediately.

2.9 Only use original parts

- Only use accessories and optional articles named or authorised by the manufacturer.
- Only use only original wear parts and replacement parts.



The manufacturer and distributor accept no liability for damages or injury resulting from the use of non-approved accessories, optional accessories, or from the use of non-original wear parts or replacement parts.

The use of non-approved accessories, optional accessories or non-genuine wear parts / replacement parts (e.g. mains cables) can have a negative effect in terms of electrical safety and EMC.

2.10 Transport

The original packaging provides optimum protection for the unit during transportation.

If required, the original packaging for the unit can be ordered.



The manufacturer and the distributor do not accept liability, even during the warranty period, for damage during transportation due to improper packaging.

- Only transport the unit in its original packaging.
- Keep the packing materials out of the reach of children.

2.11 Disposal



The unit may be contaminated. Instruct the company disposing of the waste to take the relevant safety precautions.

- Decontaminate potentially contaminated parts before disposing of them.
- Uncontaminated parts (e.g. electronics, plastic and metal parts etc.) should be disposed of in accordance with the local waste disposal regulations.
- If you have any questions about the correct disposal of parts, please contact your dental trade supplier.



An overview of the waste keys for Dürr Dental products can be found in the download area:



<http://qr.duerdental.com/P007100155>

2.12 Cybersecurity measures (security context)

The device is connected to an IT network that can be connected to the Internet. Therefore, the security context, which is a shared responsibility between the manufacturer and operator in terms of cybersecurity, must be taken into account.

The following **administrative measures** are required:

- Train employees in security and data protection.
- Appoint a data protection officer.
- Appoint a medical IT risk manager.
Use of the IEC 80001-1 standard is recommended to manage security risks for medical IT networks.

The configuration of a medical IT network requires the expertise of a network administrator. The following changes to a medical IT network can create new security risks that require a reassessment of the risk analysis:

- Changes to the IT network configuration
- Adding components (hardware or software) to the medical IT network
- Removing components from the medical IT network
- Updating or modernizing components in the medical IT network



Report cybersecurity incidents to the manufacturer at incidents@duerrdental.com.

The following **technical measures** must be taken by the operator in order to use the device safely:

- Protect IT network with firewall and antivirus software.
- Look for evidence of possible virus infection and, if applicable, check with the antivirus software and remove the virus.
- Regularly update the operating system, firewall, antivirus software and all other software applications in the IT network.
- Make sure that only trustworthy content is downloaded. Only install software and firmware updates that have been authenticated by the manufacturer.
- Segment IT network.
- Secure hardware against unauthorized access (IT network, server and connected devices).




For more information, see “Cybersecurity White Paper”.
<http://qr.duerdental.com/2110100055>



Product description

3 Overview



- 1 Combination suction unit
- 2 Mains cable with local mains plug
- 3 Hose LW 50, 1.5 m
- 4 Hose LW 50, 0.6 m
- 5 Waste water hose LW 20
- 6 Set of connection fittings
- 7 OroCup

3.1 Scope of delivery

The following items are included in the scope of delivery (possible variant-specific deviations due to country-specific requirements and/or import regulations):

Tyscor VS 4 **7188100100**

Tyscor VS 2 Plus **7188200100**

- Combination suction unit
- Mains cable
- Hose LW 50 (0.6 m)
- Hose LW 50 (1.5 m)
- Waste water hose LW 20
- Set of connection fittings
- OroCup
- Short information

3.2 Optional items

The following optional items can be used with the device:

Wall bracket 7188100015

Bracket for floor-mounted installation 7188100017

Bacteria filter 0705-991-50

Noise reduction for exhaust air . . . 0730-991-00

Flow accelerator 7560-992-00

Screed frame for flow accelerator . 7560-993-00

Surge tank 7130-991-51

Radial blower 7188100060

Rinsing unit conversion kit for Tyscor 7186100400

Upgrade hit

V/VS 2 Plus to V/VS 4 7188100075

Installation set for Tyscor V/VS

Tandem 7188100300

Tyscor Tandem Flex installation kit . 7188100302

Tyscor Tandem Flex / VS 600 installation kit (up to build year

03/2022) 7188100303

Suction flow measuring device . . . 0700-060-51

3.3 Consumables

The following materials are consumed during operation of the device and must be ordered separately:

Orotol plus (2.5 litre bottle) CDS110P6150

Orotol plus pH 7 (2.5-litre bottle) . CDS117A6150

MD 555 cleaner (2.5 litre bottle) . CCS555C6150

MD 555 cleaner organic (2.5-litre bottle) CCS556A6150

3.4 Wear parts and replacement parts

The following working parts need to be changed at regular intervals (refer to the "Maintenance" section):

Nonreturn valve (pack of 3) 7128-100-03E

Backflow preventer valve (pack of 2) 7188100082



Information about replacement parts can be found on the website portal for specialist dealers under:
www.duerrdental.net



If the mains cable of this unit is damaged it must only be replaced by an original mains cable from the manufacturer.

4 Technical data

Electrical data		7188100100	7188200100
Rated voltage	V	230, 1~	230, 1~
Mains frequency	Hz	50 / 60	50 / 60
Nominal current	A		6.2
Active power	kW		1.4
Fuses		2x T 10.0 AH / 250 V~ (IEC 60127-2)	
Type of protection			IP 21
Protection class			I

Control connection electrical data

Output:		
Voltage	V DC	24
Max. current	mA	160
Input impedance	kΩ	6.9
Hi level	V	10-30
Lo level	V	0-2.5

Electrical data for connection of the rinsing unit

Voltage	V DC	24
Max. current	mA	300

Connections

Vacuum connection (external)	mm	Ø 50
Exhaust air connection (external)	mm	Ø 50
Waste connections (DürrConnect)	mm	Ø 20

Media

Max. number of users		4	2
Max. rate of flow of fluids	l/min		20
Max. suction height	cm		80

General data

Radial blower speed (n_v) max.	At least ⁻¹		24000
Separation stage speed (n_s)	At least ⁻¹		2780 / 3170
Duty cycle	%		100
Heat generation rate	MJ/h		5.04
Dimensions (H x W x D)	cm		51 x 45 x 69
Weight	kg	24	24
Sound pressure level * approx.	dB(A)	61 / 64	61 / 64



General data

* Sound pressure level in accordance with ISO 3746

Network connection

LAN technology	Ethernet
Standard	IEEE 802.3u
Data rate	Mbit/s
Connector	RJ45
Type of connection	Auto MDI-X
Cable type	≥ CAT5

Ambient conditions during storage and transport

Temperature	°C	-10 to +60
Relative humidity	%	< 95

Ambient conditions during operation

Temperature	°C	+10 to +40
Relative humidity	%	< 70
Altitude above mean sea level	m	< 2000

Classification

Medical Device Class (MDR)	Ila
----------------------------	-----

Classification in accordance with EN ISO 10637

Classification based on separation of solids and liquids	Semi-dry suction system
Classification based on flow rate	Type 1

4.1 Characteristic curves




Fig. 1: Characteristic curves for 7188200100, measured in accordance with ISO 10637





Fig. 2: Characteristic curves for 7188100100, measured in accordance with ISO 10637

4.2 Type plate

The type plate can be found on the housing lower part.



1 Type plate


4.3 Evaluation of conformity

This device has been subjected to conformity acceptance testing in accordance with the current relevant European Union guidelines. This equipment conforms to all relevant requirements.

5 Function

When treating their patients with dental instruments, dentists, dental assistants and hygienists all face increased exposure to infectious aerosols and droplets. Studies have shown that intraoral spray mist extraction can effectively eliminate these aerosols and droplets. According to one study, a flow rate of at least 300 l/min reduces particle emissions when using a turbine to below the limit of detection (study by Koch & Graetz).

In order to achieve this volume flow at the suction cannula, a type 1 suction unit (in accordance with DIN EN ISO 10637) with an air volume flow of at least 250 NL/min is required.



- | | |
|---|--|
| 1 Network connection | 9 Orange LED – Ethernet |
| 2 Mains connection (power plug) | 10 Red LED – fault |
| 3 Connection to amalgam separator (appliance power inlet) | 11 Green LED – ready for operation |
| 4 CAN bus | 12 Blue LED – start signal |
| 5 Filter at the device suction connection | 13 Manual start button |
| 6 Intake connector | 14 SD card slot (Micro SD) |
| 7 Waste water connection | 15 Control connection |
| 8 Exhaust air connection | 16 Voltage supply for the rinsing unit |
| | 18 Touch screen |

The suction machine is a combination suction unit that is used in “wet” suction systems and consists of a separation system and a suction unit with two radial blowers. These are each driven by their own motors. The suction unit is switched on when a suction hose is removed from the hose manifold.

With the optional upgrade kit, the Tyscor VS 2 Plus can be upgraded to the Tyscor VS 4 to support up to four practitioners.

5.1 Separation system

In the separation system the aspirated fluids and the solid particles are separated from the suction air.

The waste water pump in the separation system feeds the fluid together with the finer solid particles through the waste water system connection into the central waste water network.

The speed of the separation system is continuously monitored to detect increased contamination or flooding. If the speed falls below a minimum level, the centrifugal blowers are deactivated before a fault can occur.

5.2 Radial blower

The air that has been separated from the fluids is sucked into the radial blowers. The motors in the radial blowers are regulated on a demand-driven basis by the unit electronics. Afterwards, the aspirated air is passed through the exhaust air connections and out of the device.

The antibacterial treatment of the surfaces of the radial blower with silver phosphate glass helps impede the growth of bacteria for a hygienic surface of the components.

5.3 Bacteria filter (optional)

The optional exhaust air bacteria filter reduces particles and microorganisms from the exhaust air.

For hygienic reasons, we recommend the installation of a bacteria filter in the exhaust air line.




If the device is installed in the dental practice and the exhaust air cannot be discharged to the outdoors, a bacteria filter must be installed.

5.4 LEDs

The LEDs show the status of the device:

- Operating status (ready for operation or Eco-Stop)
- Start signal of the treatment unit
- Network connection
- Fault

5.5 Touch screen



The touch screen shows the current status along with messages on the device.

Configuration is also carried out on the touchscreen.

5.6 Operating modes

Three different operating modes are available for Tyscor V/V/S 2 Plus and Tyscor V/V/S 4. These can be set up on the touch screen. If you are using monitoring software then you can also change the settings there.

The operating modes are "Balanced", "Boost" and "Eco". The different operating modes can be used to adjust the available power to the suction system (routing of lines, line lengths, layout of treatment unit etc.).

5.7 Quick Start

The Quick Start function is used to hold the unit in standby mode for a certain length of time, which ensures fast availability of the suction flow when required. Since the motor remains switched on at a very low speed when the unit is in standby mode, this option consumes slightly more energy.

5.8 Lag time


The device will continue to run for a few seconds after the suction hose has been hung back in the manifold. The lag time ensures that the remaining fluid in the suction system can be processed. Depending on the installation, the lag time can be adjusted accordingly.

5.9 Eco Stop

The Eco Stop function is used to protect the unit if it is operated inadvertently with no flow rate or with a flow rate that is too low. If the unit is operated under these conditions without the manifold signal being actuated in the mean time, the unit

will switch off automatically after a pre-defined period of time (it is possible to set this up so that it can be adjusted via the monitoring software). To switch it on, lift a suction hose up out of the hose manifold or briefly disconnect the power supply from the unit.

5.10 VistaSoft Monitor (optional)



- 1 Local practice network
- 2 Computer in the local network with server installation

The device can be connected to the VistaSoft Monitor monitoring software via the local practice network.

The software displays the device's messages, current status and upcoming maintenance. The information can be viewed on the computer in the practice or via the cloud.

5.11 Grouping of multiple devices


By combining multiple devices together in this way, the available volume flow can be increased, which in turn means that the maximum number of therapists can use the system.

Up to three Tyscor suction machines can be connected together with the *installation set for Tyscor V/VS Tandem*.

The *installation set for Tyscor Tandem Flex* is required to combine a Tyscor suction unit with a side channel suction unit.

See also "3.2 Optional items".

Tyscor Tandem: Combination of several Tyscor suction units




Up to three Tyscor suction units can be operated together as a group. Here, one device takes on the role of the controller (main unit) and controls up to two additional devices (auxiliary units).



When grouped together, mutual controlling and exchange of information takes place via the CAN bus. Only one device is permitted to be configured as the main unit; additional devices must be configured as auxiliary units.


Tyscor Tandem Flex: Combination of one Tyscor with a side channel suction unit



The Tyscor suction unit can be operated in combination with a side channel suction unit. The Tyscor monitors the vacuum in the suction line and switches itself on and off automatically as required. It thus serves to extend the side channel suction unit.

In this constellation, only the "Balanced" and "Boost" operating modes are available.

5.12 OroCup



The OroCup care system is a closed dosing system for simple preparation and suctioning of disinfectants and special cleaners.

Use OroCup to clean and disinfect the suction system with all of its components and the spittoon.



6 Requirements

The unit can be installed on the same level as the surgery room or in a floor below.



Further information can be found in our suction planning information leaflet.
Order number 9000-617-03/..

6.1 Installation/setup room

The room chosen for set up must fulfil the following requirements:

- Closed, dry, well-ventilated room
- Should not be a room made for another purpose (e. g. boiler room or wet cell)
- When installing in a cabinet the inlet and outlet ventilation slots must be present; minimum free cross-section at least 120 cm².
- Forced ventilation (fan) must be provided if there is a risk that the recommended room air temperature could be exceeded. The air flow performance must be at least 2 m³/min.
- Do not cover cooling slots or openings with housing installations; ensure sufficient clearance to the openings to permit sufficient cooling.
- Mains cable plug connections must be freely accessible so they can be quickly disconnected if there is any danger.

6.2 Setup options

The following options for setting up the unit are available:

- On the floor with a console for floor-mounted installation
- Wall installation using a Dürr Dental wall mounting
- In a ventilated cabinet

6.3 Pipe materials

Only use vacuum-sealed HT-waste pipes manufactured from the following materials:

- Polypropylene (PP),
- Chlorinated polyvinyl chloride (PVC-C),
- Plasticizer-free polyvinyl chloride (PVC-U),
- Polyethylene (PE).

The following materials must not be used:

- Acrylonitrile-butadiene-styrene (ABS),
- Styrene copolymer blends (e.g. SAN + PVC).

6.4 Hose materials

For waste connections and suction lines only use the following hose types:

- Flexible spiral hoses made of PVC with integrated spiral or equivalent hoses
- Hoses that are resistant to dental disinfectants and chemicals



Plastic hoses will display signs of ageing over time. Therefore, they should be inspected regularly and replaced as necessary.

The following types of hoses must not be used:

- Rubber hoses
- Hoses made completely of PVC
- Hoses that are not sufficiently flexible

6.5 Information about electrical connections

- Ensure that the electrical connections to the mains power supply are established in accordance with current valid national and local regulations and standards governing the installation of low voltage units in medical facilities.
- Observe the current consumption of the devices that are to be connected.

6.6 Information about connecting cables


Mains supply cable

Only use the supplied mains cable to connect the device.

Control cable

Installation type	Line layout (minimum requirements)
Fixed installation	– Shielded sheathed cable (e.g. (N)YM (St)-J)

Installation type	Line layout (minimum requirements)
Flexible	<ul style="list-style-type: none"> – PVC data cable with shielded cable sheathing, as used for telecommunications and IT processing systems (e.g. type LiYCY) <p>or</p> <ul style="list-style-type: none"> – Lightweight PVC control cable with shielded cable sheathing

 Connect the shielding of the cables in accordance with the regulations.

7 System components


The system components listed below are required or recommended for various procedures or for installation.

7.1 Rinsing unit

It is recommended that the suction system is equipped with a rinsing unit, e.g. in the treatment unit. The rinsing unit provides a small amount of water during aspiration. This dilutes the aspirated fluids (blood, saliva, rinsing water etc.), which can then be transported more effectively.

7.2 Exhaust air filter / bacteria filter

For hygienic reasons, we recommend the installation of a bacteria filter in the exhaust air line. If the unit is installed in the surgery and the exhaust air cannot be discharged to the outdoors, it is essential to install a bacteria filter. Depending on the type and condition of the bacteria filter, it will need to be replaced every 1-2 years at the latest.

 The separation integrated in the system does not retain bacteria; this is why we recommend installing a suitable filter in the exhaust air system.

7.3 Noise reduction

If the noise level from the exhaust air vent or the flow noise generated is too high, noise reduction can be installed in the exhaust air line.

7.4 Surge tank

If the suction unit is combined with an amalgam separator, this requires the installation of a surge tank. The surge tank reduces pressure peaks caused by the waste water pump of the suction unit and acts as a buffer against temporary rises in the volume of water.

The surge tank can also be used if the waste water is fed directly into the building waste water system. In this case the waste water from the suction unit is diverted to the building drainage system under zero pressure.

7.5 Flow accelerator

In order to keep the suction system free of deposits, a flow accelerator can be fitted in conjunction with a spittoon valve. When using a bowl rinse system, water will collect before the flow

accelerator. The next time suction takes place using the large cannula, the collected fluid is transported in surges and at high speed to the suction system. This ensures automatic cleaning of the suction pipes.


8 Installation



The actual connection can vary depending on the chosen installation option. The connection shown is only an example.

8.1 Installation and routing of hoses and pipes

- Establish connections between the pipe system and the unit using the flexible hoses supplied. This will prevent vibrations from being transmitted to the pipe system.
- The connection between the pipe line and unit suction connection should be kept as short as possible and straight, without bends.
- Install the drain hoses with a downward gradient so that the waste water can drain off.
- Waste water connections must be implemented in accordance with applicable local and national regulations.



- 1 Curved adaptor bush
- 2 Straight adaptor bush
- 3 Sealing sleeve
- 4 Union nut
- 5 Suction hose Ø 50 mm (internal)
- 6 Hose clamp Ø 55 mm
- 7 Exhaust air hose Ø 50 mm (internal)
- 8 Hose clip Ø 28 mm
- 9 Hose sleeve
- 10 Waste water hose Ø 20 mm (internal)
- 11 Securing ring
- 12 Hose sleeve Ø 20 mm
- 13 O-ring 20 x 2.0
- 14 Connector Ø 36 mm (external)
- 15 O-ring 30 x 2

8.2 Installing the group of devices

Various mounting kits are available for installing several devices in a device network (see also "3.2 Optional items"):

- *Installation set for Tyscor V/VS Tandem* for the combination of up to three Tyscor suction machines
- *Installation set for Tyscor Tandem Flex* mounting kit for combining a suction unit with a Tyscor side channel suction machine

The installation sets contain various parts (pipes, fittings, hose clamps, etc.) and corresponding installation information.




Setup and installation must be carried out in accordance with the requirements of Dürr Dental. Corresponding installation information is enclosed with the installation set but can also be downloaded from the Dürr Dental homepage.

The *installation set for Tyscor Tandem Flex / VS 600 (up to production date 03/2022)* is also required for connection to a VS 600 side channel suction unit.



All suction units in the group must be on the same level, and the supply manifold must be at least 5 cm lower than the suction inputs.



- 1 Rinsing unit
- 2 Water connection
- 3 Voltage supply for the rinsing unit
- 4 Rinsing hose
- 5 Overflow connection piece

8.3 Fitting the rinsing unit

1. Securely screw the rinsing unit onto the retaining plate.
2. Securely screw the rinsing unit and the retaining plate onto the housing upper part.
3. Connect the hose for the water supply to the rinsing unit and to the water supply.
4. Push the rinse hose onto the rinse connection of the separation unit.
5. Shorten the rinse hose to the required length and connect it to the rinsing unit.
6. Connect the waste water hose to the overflow connection piece and to the waste outlet.
7. Connect the voltage supply of the rinsing unit at the suction unit.


8.4 Network connection

Purpose of the network connection

The network connection is used to exchange information or control signals between the unit and a software installed on a computer, in order to, e. g.:

- Display parameters
- Select operating modes
- Indicate messages and error situations
- Change unit settings
- Activate test functions
- Transmit data for archiving
- Provide documents concerning the units

1. Plug in the network cable at the network connection on the unit (optional when using monitoring software).
2. Plug in the network cable at the network socket.




1 Network cable

8.5 CAN bus

Purpose of the CAN bus connection

In a grouping of up to three devices, the CAN bus is used to enable mutual controlling and information exchange between the devices.

1. Connect the network cable between the devices at the CAN bus.



- 1 Main unit
- 2 Auxiliary unit 1
- 3 Auxiliary unit 2


8.6 Electrical connections




WARNING

Electric shock


- › The device may only be connected to a supply system with a earthed power outlet.

1. Fasten the plug socket to the control line and connect to the device.
 In a grouping of devices connect the control cable to the main unit.
2. Connect the mains cable to the unit and to the power outlet.



- 1 Mains cable with socket and country-specific mains plug
- 2 Control line

8.7 PCB (main board) electrical connections



- | | |
|-----|--|
| X1 | Supply voltage 230 V |
| X2 | Control voltage output, 24 V DC, control voltage input |
| X3 | Voltage supply for the rinsing unit, 24 V DC |
| X5 | SD card holder (for Micro SD) |
| X6 | CAN bus |
| X7 | Service interface |
| X8 | Connection for the display |
| X10 | Supply voltage for suction motor 1 |
| X11 | Supply voltage for suction motor 2 |
| X12 | Separation motor supply voltage (VS only) |
| X13 | Separation motor RPM monitor (VS only) |
| X14 | Fan connection for housing fan |
| X15 | Motor control fan connection 1 |
| X16 | Motor control fan connection 2 |
| X17 | Connection for cooling air fan (optional) |
| X18 | Jumper (V = closed, VS = open) |
| X20 | Network connection |
| F1 | Mains fuse T10.0AH |
| F2 | Mains fuse T10.0AH |
| H1 | Blue LED – start signal |
| H2 | Red LED – fault |
| H3 | Green LED – voltage supply +24 V |
| H4 | Green LED – voltage supply +3.3 V |
| H5 | Green LED – voltage supply +5 V |

H6	Green LED – HV voltage supply +15 V
H7	Green LED – HV voltage supply +5V
H8	Green LED – HV voltage supply +3.3 V
H9	Green LED – temperature indicator for radial blower 1, temperature OK
H10	Yellow LED – line to radial blower 1, short circuit or open circuit
H11	Red LED – temperature indicator for radial blower 1, temperature too high
H12	Green LED – temperature indicator for radial blower 2, temperature OK
H13	Yellow LED – line to radial blower 2, short circuit or open circuit
H14	Red LED – temperature indicator for radial blower 2, temperature too high
H21	Green LED – ready for operation
H22	Orange LED – network
S1	Start button
B10	Vacuum sensor
E1	Battery (CR2032 button cell)

8.8 Testing suction power after installation

Measure negative pressure

The following procedure is recommended in order to measure the maximum negative pressure of the suction system so that it can be compared with the performance curve:

1. Select the operating mode "Boost"
2. At a treatment unit lift off a suction hose and connect a pressure gauge (relative pressure). The connection adapter to the pressure gauge must be leak tight.
3. Read off the displayed negative pressure and compare it against the performance curve.

Measure volume flow

To ensure that all particles and liquids sucked in are reliably removed, 300 - 350 l/min must be extracted at the suction cannula.

The volume flow can be measured with a suction flow meter, see "3.2 Optional items".

See also the installation and operating instructions for the measuring device.

1. Check the volume flow at the large suction hose on the treatment unit.
Simultaneously lift the large suction hose depending on the maximum number of users on the treatment units.

If the volume flow is lower, the following causes may be present in the suction unit:

- Diameter of the pipe system is too small or unfavourable pipe construction
- Hose cross-sections are too small (instrument hose or hoses in the treatment unit)
- Dirty coarse/fine filter
- Hoses clogged
- Suction unit underpowered or defective

9 Commissioning



In many countries technical medical products and electrical devices are subject to regular checks at set intervals. The owner must be instructed accordingly.



NOTICE

Interference caused by larger particles such as pieces of tooth or fillings

- › Do not operate the unit without a coarse filter

1. Check that the coarse filters are installed in the suction system (e.g. in the spittoon).
2. Turn on the unit power switch or the main surgery switch.
3. Carry out a function check of the device.
4. Check all connections for leak tightness.
5. Carry out an electrical safety check in accordance with applicable regulations (e.g. regulations concerning set up, operation and application of medical devices) and record the results as appropriate, e.g. in the technical log book.
6. Carry out and document the instruction and handover for the unit.



A sample handover report is included in the attachment.


9.1 Setup wizard

A setup wizard is displayed on the touch screen when the device is started for the first time or after a firmware update. To ensure trouble-free operation, please follow the instructions displayed in the setup wizard.

9.2 Settings on the touch screen

In the Settings screen you can query information and adjust settings.

To access the individual menu items, tap the relevant button. If there are multiple pages, the arrow buttons are used to scroll through the pages.



Changing to the Administrator user level:

In the display screen, use the **Settings** button to switch to the menu screen. Use the **Access Level** button to switch to the selection for **User**, **Administrator** and **Service Technician [PIN]**. Then tap **Administrator** and confirm. Afterwards use the **Main Menu** button to go back.

Selecting the operating mode

Depending on the installation, the operating mode (suction power) of the device can be adjusted. Before changing the operating type, the suction power should be checked at the can-nulas.

The operating mode can be adjusted in the display screen by tapping the buttons **Eco**, **Bal-anced** or **Boost**. The selected operating mode is highlighted with a different colour.

Adjust Tyscor Tandem group

If several Tyscor suction units are operated in a group, one of the suction units must be defined as the **main unit**. Each additional suction unit is then defined as **auxiliary unit 1** or **auxiliary unit 2**.



On delivery, each device is preconfigured as a **main unit**.

Requirements:

- ✓ The suction units must be connected to each other via CAN bus.
- ✓ You need to be logged in as an **Administrator** on the touchscreen.

Setting up the group settings:

1. Use the **Settings** button to switch from the display screen to the menu screen.
2. Use the **Parameters** button to switch to **Cluster Settings**.
3. In the selection list, select whether the device is a **main unit** or an **auxiliary unit**.

Control strategy in group operation

If required, you can change the control strategy of the device or device group.

Setting up the control strategy:

1. Use the **Settings** button to switch from the display screen to the menu screen.
2. Use the **Parameters** button to go to **Control strategy options**.
3. In the selection list, select the desired control strategy of **Demand-based Mode** or **Parallel Mode**



In **Parallel Mode**, increased power consumption and increased wear are to be expected.

Adjust Tyscor Tandem Flex group

In order to be able to operate the Tyscor suction unit in combination with a side channel suction unit, it is necessary to select the control strategy **Tyscor Tandem Flex** on the Tyscor suction unit.

Requirements:

- ✓ You need to be logged in as an **Administrator** on the touchscreen.

Setting up the control strategy:


1. Use the **Settings** button to switch from the display screen to the menu screen.
2. Use the **Parameters** button to go to **Control strategy options**.
3. Select the control strategy **Tyscor Tandem Flex** and save.



In this operating mode it is only possible to select the modes **Balanced** and **Boost**.

Setting the afterrun delay time

Depending on the installation, the lag time of the device can be adjusted. The lag time can only be adjusted from the user level **Administrator** or higher.


 When operating devices as a group, always adjust the lag time on the main unit.

Adjusting the lag time:

To adjust the lag time, switch from the display screen to the menu screen via the **Settings** button. Then use the **Parameters** and **Lag Time** buttons to adjust the time. Tap the **Save** button to save the time.

Quick Start

Depending on the installation, the startup time of the device can be adjusted. The startup time can be adjusted from the user level **Administrator** or higher.


 When operating devices as a group, set up the Quick Start function on the main unit.

Activating Quick Start:

To activate the Quick Start function, switch from the display screen to the menu screen by tapping the **Settings** button. Then switch to the submenu via the buttons **Parameters** and **Start/Stop**. The function can now be enabled or disabled via the **Quick Start** slider. If the function is active then the slider is shown in blue.

Eco Stop

The Eco Stop function can be set to either a short time (enabled) or long time (disabled). This setting can be adjusted from the user level **Administrator** or higher.

 When operating devices as a group, set up the Eco Stop function on the main unit.

Enabling the Eco Stop:

To activate the Eco Stop function, switch from the display screen to the menu screen by tapping the **Settings** button. Then switch to the submenu via the buttons **Parameters** and **Start/Stop**. The function can now be enabled or disabled via the **Eco Stop** slider. If the function is active then the slider is shown in blue.

Exhaust air filter

If an exhaust air filter is used, this will need to be set up in the controller. This is the only way to ensure that a corresponding maintenance message is displayed at the right time.

Settings for the exhaust air filter:

To activate the maintenance counter for the exhaust air filter, switch from the display screen to the menu screen by tapping the **Settings** button. Then use the **Parameters** and **Exhaust Air Filter** buttons to switch to the selection list. In the list, select whether or not an exhaust air filter is fitted.

Emergency mode

In emergency mode the device or group of devices can continue to be operated even if one of the device components is defective. Once the device component has been repaired, emergency mode will need to be reset. This can be done from user level **Administrator** or higher.

Resetting emergency mode:

To reset emergency mode, switch from the display screen to the menu screen with the **Settings** button. Then use the **Service** and **Emergency Mode** buttons to switch to the submenu. Use the button to reset emergency mode.

Device components


In the **Device Components** menu you can deactivate one or possibly both radial blowers for maintenance purposes. The menu is available beginning in the **Administrator** user level.

 Both of the radial blowers must be activated in normal operation.

Deactivating/activating radial blower(s):

To deactivate/activate radial blowers, switch from the display screen to the menu screen with the **Settings** button. Then use the **Service** and **Device Components** buttons to switch to the submenu. The desired motor can be deactivated/activated via the sliders for the top motor and bottom motor. If the motor is active then the slider is shown in blue.

During operation, the status of the motors is displayed in colour in the **Device Components** submenu. A green text shows that the motor is running. A red text shows that the motor has a fault.

 Status changes are only displayed after calling up the **Device Components** submenu.

9.3 Monitoring the unit via the network

If the device is connected to the VistaSoft Monitor monitoring software via the local practice network, the software displays the device's messages, current status and upcoming maintenance.

The information can be viewed on the computer in the practice or via the cloud.

For further information, see VistaSoft Monitor handbook:



<http://qr.duerdental.com/2110300001>

The following requirements must be met in order to monitor the unit on the computer:

- Unit connected to the network
- Current monitoring software installed on the computer

Combining devices safely

- The overall safety of the unit and its main performance features are independent of the network. The device is designed for operation independent of a network. However, some of the functions are not available in this case.
- Incorrect manual configuration can lead to significant network problems. The expert knowledge of a network administrator is required for configuration.
- The data connection utilises part of the bandwidth of the network. Interactions with other medical devices cannot be completely ruled out. Apply the IEC 80001-1 standard for risk assessment.
- The device is not suitable for direct connection to the public Internet.
- Protect the network with antivirus and firewall tools.
- Only allow authorised users to access the devices.
- Only allow authorised users to access the device network.
- Block access to the device network via the guest WiFi.
- Only download firmware updates from recognised, secure sources, e.g. directly from the manufacturer's website.

Network configuration

Various options are available for network configuration:


- ✓ Automatic configuration via DHCP (recommended).
 - ✓ Automatic configuration via Auto-IP for direct connection of unit and computer.
 - ✓ Manual configuration.
1. Configure the network settings of the unit using the software or, if available, the touch screen.
 2. Check the firewall and release the ports, if applicable.

Network protocols and ports


Port	Purpose	Service
45123 UDP, 45124 UDP	Unit recognition and configuration	
1900 UDP	Service indicator	SSDP / UPnP
502 TCP	Device data	
514 ¹⁾ UDP	Event log data	Syslog
22 TCP, 23 TCP	Diagnosis	Telnet, SSH
123 UDP	Time	NTP







- ¹⁾ The port may vary depending on the configuration.

Usage

 When using prophylaxis powders, water-soluble Lunos Prophylaxis Powders are recommended in order to protect the suction systems (Dürr Dental).

10 LEDs




- | | | | |
|---|----------------------------|---|--|
| 1 | Network |  | <p>ORANGE display is illuminated</p> <p>The machine is connected to the network.</p> |
| 2 | Error |  | <p>Red display flashes</p> <ul style="list-style-type: none"> – Fault/failure of a radial blower |
| | |  | <p>RED LED illuminates</p> <ul style="list-style-type: none"> – Failure of both radial blowers or – Failure of separation (VS only) |
| 3 | Operation |  | <p>Ready for operation</p> <p>Green LED lights up</p> |
| | |  | <p>EcoStop</p> <p>Green LED flashes</p> <p>The unit has been switched off by the EcoStop function. To switch it on, lift a suction hose up out of the hose manifold or briefly disconnect the power supply from the unit.</p> |
| 4 | Hose manifold start signal |  | <p>BLUE LED illuminates</p> <p>Manifold signal active and machine running.</p> |

11 Operating the touch screen

The selected operating mode and the current negative pressure are displayed on the Main Menu of the main unit. The operating mode can be changed by tapping one of the operating modes.





The **Next** button can be used to scroll to the next page, where the power consumption and temperature of the electronics are displayed.

The power consumption and temperature of the electronics are displayed on the Main Menu of the auxiliary unit.




- 1 Context range (e.g. date, time)
- 2 Content area
- 3 Navigation buttons

The buttons at the bottom edge of the screen can be used to switch between the different menus.

-  Switches to the menu **Settings**
-  Switches between the two start screens (only when a device is configured as the "main unit")
-  Switches to the start screen
-  Switches to the next higher menu level





11.1 Querying messages

Any incidents that occur in the unit are shown in the Messages screen. The messages can be confirmed by tapping the navigation button. The messages are then displayed in the context area of the Home screen until the incident has been remedied.



- 1 Scrolling between messages
- 2 Detailed information about the message
- 3 Navigation buttons

Severity of the messages:

-  Information
-  Note
Important information about the device
-  Notice
Operation of the device restricted
-  Fault
Operation of the device interrupted

12 Disinfection and cleaning

The following tasks are required for the disinfection and cleaning of the suction system:

"12.1 Suctioning water"	After every treatment
"12.2 Disinfection"	Daily in the evening after the end of treatment, With increased workloads, before the midday break and after the end of treatment
"12.3 Cleaning"	Allow to act for at least 2 x per week before the start of treatment or during the midday break Alternatively, 5 x per week, with short reaction times




NOTICE

Device malfunctions or damage due to use of incorrect media

Guarantee claims may become invalid as a result.

- › Do not use any foaming agents such as household cleaning agents or instrument disinfectants.
- › Do not use abrasive cleaners.
- › Do not use agents containing chlorine.
- › Do not use any solvents like acetone.

The following should always be used:



- For disinfection and cleaning:
 - Orotol plus
 - Orotol plus pH 7
- For cleaning:
 - MD 555 cleaner
 - MD 555 cleaner organic


Only these products have been tested by Dürr Dental.




Comply with the instructions for use of the cleaning agent and disinfectant.

12.1 Suctioning water


After every treatment:

1. Suction up cold water (min. 0.5 litres) with the large and small suction hoses. Do this even if only the small suction hose was actually used during treatment.
-  Suction through the large suction hose causes a large amount of air to be drawn up, thereby considerably increasing the cleaning effect.





12.2 Disinfection


Disinfect and clean the suction system every evening after the end of treatment.

-  With increased workloads, clean and disinfect twice per day, e.g. before the mid-day break and after the end of treatment.

The following is required for the disinfection and cleaning of the suction system:

-  Non-foaming disinfectant/cleaning agent that is compatible with the materials
 - Orotol plus
 - Orotol plus pH 7

-  care system
 - e.g. OroCup

-  Disinfection and cleaning are described below with OroCup and Orotol plus.



Wear protective gloves.




Wear protective goggles.

1. To pre-clean, suction up 2 litres of water with the care system.




Using the solution in the care system:

1. Pour 2 sealing caps (40 ml) of Orotol plus into the OroCup.



2. Fill the OroCup with 2 litres of cold water.



3. Close the lid of the care system.



4. Mix the solution.



Using the solution:

1. Attach suction hoses to the care system and suction up 1 litre of solution.



2. Pour the remaining solution into the spittoon.



3. Allow to work for a minimum of one hour or leave overnight.

4. When placing the system back into operation, suction up 2 litres of water.



The instructions are also available as a video:



12.3 Cleaning

Clean the suction system regularly:

- Allow to act for at least 2 x per week in the morning before the start of treatment or during the midday break, reaction time 30–120 minutes
- Alternatively: 5 x per week before the start of treatment, reaction time 5–15 minutes

The following is required for cleaning the suction system:



Special non-foaming detergent for suction systems that is compatible with the materials

- MD 555 cleaner
- MD 555 cleaner organic



care system
– e.g. OroCup



Cleaning is described below with OroCup and MD 555 cleaner.



Wear protective gloves.




Wear protective goggles.


1. To pre-clean, suction up 2 litres of water with the care system.



2. Use the solution in the care system.
Pour 5 caps (100 ml) of MD 555 cleaner into the care system.



3. Fill the OroCup with 2 litres of cold water.



4. Close the lid of the care system.



5. Mix the solution.



6. Attach suction hoses to the care system and suction up one litre of solution.



7. Pour the remaining solution into the spittoon.



8. Allow the solution time to act.
For 2 x per week: 30–120 minutes
For 5 x per week: 5–15 minutes

9. Suction up 2 litres of water after the reaction time has elapsed.



The instructions are also available as a video:



13 Maintenance



All maintenance work must be performed by a qualified expert or by one of our Service Technicians.

Maintenance interval	User maintenance work
Every 3 months	› Check the filter at the device suction connection and clean if necessary.
Maintenance interval	Maintenance work by service technician
Annually	› Replace the waste water valve. › Check the nonreturn valves upstream of the radial blowers for correct operation and replace as required.
Every 1-2 years or when prompted if the maintenance counter is activated	› Replace the exhaust air filter (where fitted).



Note maintenance information on the touchscreen and, if applicable, in VistaSoft Monitor.



WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).



Prior to working on the unit or in case of danger, disconnect it from the mains.

13.1 Wear parts and replacement parts

The following wear parts must be replaced at regular intervals:

Nonreturn valve (pack of 3)	7128-100-03E
Backflow preventer valve (pack of 2)	7188100082



Information about replacement parts can be found on the website portal for specialist dealers under:
www.duerrdental.net

13.2 Clean the filter at the device suction connection.

Clean the filter at the suction connection every 3 months.




WARNING


Infection due to contaminated unit

- Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).

1. Loosen the filter at the suction connection by twisting it half a turn anti-clockwise.



2. Clean the filter.




3. Push the filter with the opening facing upwards into the suction connection at the separation housing.
4. Screw in the protective filter clockwise as far as it will go.

13.3 Confirm maintenance

The maintenance interval can be reset on the touchscreen if maintenance has been carried out.

Requirements:

- ✓ Logged in as an Administrator on the touchscreen

1. Touch the touchscreen .
2. Touch **Maintenance**.
3. Select the maintenance task that was performed.
4. Touch **Maintenance Performed**.

13.4 Intensive cleaning

Intensive cleaning can be carried out during the self-clean process of suction systems that are heavily contaminated with biofilms or if the suction power is greatly reduced.

Carry out this intensive cleaning at least twice a week.

Exposure time: 1 – 2 hours

The following is required for intensive cleaning of the suction system:



Special non-foaming detergent for suction systems that is compatible with the materials

- MD 555 cleaner
- MD 555 cleaner organic



care system

- e.g. OroCup



Cleaning is described below with OroCup and MD 555 cleaner.



Wear protective gloves.




Wear protective goggles.

1. To pre-clean, suction up 2 litres of water with the care system.



2. Use the solution in the care system. Pour 10 caps (200 ml) of MD 555 cleaner into the care system.



3. Fill the OroCup with 2 litres of cold water.



4. Close the lid of the care system.



5. Mix the solution.



6. Attach suction hoses to the care system and suction up one litre of solution.



7. Pour the remaining solution into the spittoon.



8. Allow the solution to act for 1 – 2 hours.

9. Suction up 2 litres of water after the reaction time has elapsed.



? Troubleshooting

14 Tips for operators and service technicians



Any repairs exceeding routine maintenance may only be carried out by qualified personnel or our service.



WARNING

Infection due to contaminated unit

- › Clean and disinfect the suction before working on the unit.
- › Wear protective equipment when working (e. g. impermeable gloves, protective goggles and mouth and nose protection).



Prior to working on the unit or in case of danger, disconnect it from the mains.



WARNING

Electric shock due to capacitor discharge

- › Wait for the discharge time.
- › Watch for the LEDs going out.

14.1 General faults

Error	Possible cause	Remedy
Device does not start	No mains voltage	<ul style="list-style-type: none"> › Check the mains supply voltage. * › Check the fuses and replace if necessary. *
	Undervoltage	<ul style="list-style-type: none"> › Measure the supply voltage; call an electrician if necessary. *
	Control electronics defective	<ul style="list-style-type: none"> › Replace the electronics. *
The unit has been switched off, the green LED is flashing	The EcoStop function has been activated	<ul style="list-style-type: none"> › Lift up a suction hose from the hose manifold. › Briefly disconnect the power supply from the unit. › Check the switch in the hose manifold and replace as required. *


Error	Possible cause	Remedy
Water leaking from the exhaust air connection	Membrane valve defective	› Check the membrane valve at the waste water connection and if necessary clean or replace.
	Foam in turbine due to use of incorrect cleaning and disinfectant agents	› Use non-foaming cleaning and disinfectant agents.
	Build-up of condensate in the exhaust air line	› Check the pipe system; avoid over-cooling. *
Suction performance too low	Filter at the suction connection clogged	› Clean the filter at the device suction connection.
	Leak in the suction pipe	› Check and, if necessary, establish leak-tightness of the suction pipe and connections. *
	Soiling in the suction pipe	› Check the pipe system and clean it as required. *
	Soiling in the separation stage	› Check the separation stage and clean it as required. *
	Membrane valve defective	› Check the membrane valve at the waste water connection and if necessary clean or replace.
	One radial blower defective	› Replace the radial blower. *
No suction power	Radial blower defective	› Activate emergency mode by confirming the alarm message. › Replace radial blower. *
	Controller defective	› Replace the controller. *
	Separation system defective	› Check the separation system and clean or replace it as required. *
	Main unit defective in group	› Take the defective main unit out of operation and reconfigure an auxiliary unit as the main unit.
Water not being pumped away	Separation system defective	› Replace separation system. *
Nothing shown on the display	Cable defective or not connected	› Check the cable connection and replace the cable if necessary. *
	Display defective	› Replace the display. *

Error	Possible cause	Remedy
Display does not respond to touches	Display not calibrated	<ul style="list-style-type: none"> › Calibrate the display. – Disconnect and reconnect the mains plug. Within 8 seconds of the device being started, start the calibration process with a long press (5 s) on the display. Follow the instructions on the display.
	Display defective	› Replace the display. *
* Only to be done by service technicians.		

14.2 Fault with Tyscor Tandem Flex

Error	Possible cause	Remedy
No suction power No error message on Tyscor touchscreen. Error message on VistaSoft Monitor may appear if side channel suction unit is connected to Connect Box.	Side channel suction unit defective	› Repair or replace side channel suction unit.* › Change operation to Tyscor, see "Emergency mode: Change to Tyscor".
Low suction performance Error message on the touchscreen of the Tyscor and in VistaSoft Monitor (if available).	Side channel suction unit does not generate sufficient negative pressure, but the motor protection switch in the control box does not trip. Tyscor suction unit defective.	› Repair or replace side channel suction unit.* › Change operation to Tyscor, see "Emergency mode: Change to Tyscor". › Repair or replace Tyscor suction unit.*

* Only to be done by service technicians.


 If the suction power is greatly reduced but there is no error message on the Tyscor suction unit, remove the suction hose from the hose holder for longer than 20 seconds. In this way, the control unit recognizes the failure of the side channel suction unit and the error message is triggered.

Emergency mode: Change to Tyscor

The operation of Tyscor Tandem Flex is set on the control line switch:




0 Tyscor Tandem Flex deactivated.
The start signal is sent to the Tyscor suction unit.



I Tyscor Tandem Flex activated.
The start signal is sent to both suction units. The Tyscor monitors the vacuum in the suction line and switches on automatically if necessary.

If the motor protection in the control box of the side channel suction unit is tripped, Tyscor no longer receives a signal and so no longer switches on.

You can then switch to emergency mode in which only the Tyscor suction unit is running. Only the suction power of the Tyscor suction unit is then available.

1. Check the control box of the side channel suction unit to see if the motor protection switch has tripped. If not, switch off the control box.
2. Close the ball valve to the side channel suction unit.
3. Set the switch in the control line to "0" .
4. Confirm error message at the Tyscor.
The control strategy of the Tyscor suction unit is switched to the standard mode.
5. If no error message is displayed on the Tyscor, manually switch to standard mode on the touchscreen under **Settings > Parameters > Control Strategy Options**.

14.3 Error messages



Error messages are displayed on the touch screen. If there is a network connection, the messages can be forwarded to the monitoring software. If the device is not connected to the network, the messages can be read via a terminal client (e. g. PuTTY).

Error	Possible cause	Remedy
Maintenance required. Clean filter at device suction connection.	Service due	<ul style="list-style-type: none"> › Clean the filter at the device suction connection. › Reset the maintenance in the menu.
Maintenance required. Check waste valve.	Service due	<ul style="list-style-type: none"> › Check the waste valve. * › Reset the maintenance in the menu. *
Maintenance required. Replace exhaust air filter.	Service due	<ul style="list-style-type: none"> › Replace the exhaust air filter * › Reset the maintenance in the menu. *
Severe fault detected. Confirm to enter emergency mode.	A radial blower is defective	<ul style="list-style-type: none"> › Confirm to enter emergency mode. › Immediately repair or replace the radial blower. * › Reset emergency mode in the menu. *
Severe fault detected. Call service.	A radial blower is defective. Emergency mode unavailable.	<ul style="list-style-type: none"> › Immediately repair or replace the radial blower. *
This device is running in emergency mode. Call service.	Device is running in emergency mode	<ul style="list-style-type: none"> › Immediately repair or replace the radial blower. * › Reset emergency mode in the menu. *
An Aux device has a severe failure. Confirm at affected device.	A radial blower on an auxiliary unit is defective	<ul style="list-style-type: none"> › Confirm to enter emergency mode. › Immediately repair or replace the radial blower. * › Reset emergency mode in the menu. *
An Aux device is running in emergency mode. Call service.	An auxiliary unit is running in emergency mode	<ul style="list-style-type: none"> › Immediately repair or replace the radial blower. * › Reset emergency mode in the menu. *
Tyscor Tandem Flex operation: Side channel blower defective	Side channel fan defective	<ul style="list-style-type: none"> › Deactivate side channel fan on the control box. › Close the pipes to the side channel fan at the ball valve. › Switch the deposit signal to the Tyscor suction unit. › Confirm error message on the display.
Driver overcurrent TRIP	Radial blower motor defective	<ul style="list-style-type: none"> › Replace the radial blower. *

Error	Possible cause	Remedy
Speed of Sepa is low	Motor defective	› Replace the separation stage. *
	Hall sensor PCB (main board) defective	› Replace the Hall sensor PCB (main board), check the magnets in the Sepa fan. *
	Centrifuge soiled or damaged	› Check the centrifuge and clean or replace if necessary. *
Vacuum motor overheated	Radial blower motor defective	› Replace the radial blower. *
DC bus overvoltage	Control error	› Replace the electronics. *
DC bus undervoltage	Mains power supply fault	› Check the mains connection and supply voltage. *
	Machine was disconnected from the mains while running	› No action required.
	Control error	› Replace the electronics. *
No Ready Signal from vacuum machine	Control error	› Replace the electronics. *
Internal board communication disturbed	Failed firmware update	› Perform/repeat the firmware update. *
	Control error	› Replace electronics. *
Unexpected re-initialization	Firmware error	› Perform/repeat the firmware update. *
Short circuit to earth	Radial blower motor defective	› Replace radial blower. *
Vacuum motor sensor shorted	Radial blower motor defective	› Replace radial blower. *
Vacuum motor sensor open circuit	Radial blower motor defective	› Replace radial blower. *
	Motor cable not correctly connected to the control board	› Check the plug connection. *
Firmware mismatch	Different firmware versions on the two processors after a firmware update	› Perform/repeat the firmware update. *
Speed Feedback Failure	Motor speed detection defective	› Replace radial blower. *
MC Lib failure	Control error	› Replace the electronics. *
CPU overheated	Insufficient ventilation or poor set-up conditions	› Check the setup conditions, ensure adequate ventilation.
	Fan in the foam housing soiled	› Clean the fan and ventilation slots for supply and exhaust air. *
	Fan in foam housing defective	› Replace the fan. *
	Control electronics defective	› Replace electronics. *

Error	Possible cause	Remedy
Power Pack overheated	Insufficient ventilation or poor set-up conditions	› Check the setup conditions, ensure adequate ventilation.
	Fan on the electronics housing soiled	› Remove the cover on the electronics housing, clean the fan and heat sink. *
	Fan on electronics housing defective	› Replace the fan. *
	Control electronics defective	› Replace electronics. *
Eco Stop. Switch start signal off and on again to restart	Unit has been switched on unintentionally for too long.	› Check whether all suction hoses are correctly hung up.
		› Briefly disconnect the unit from the mains.
		› Check whether a permanent start signal is present at the unit. *
		› Check the control cable. *
All devices failed or deactivated	Radial blower defective	› Replace the radial blower. *
	Radial blower deactivated	› Inform a Service Technician.
More than one main control unit detected. Wrong configuration!	Configuration error. More than one main unit found	› Make sure that only one device is configured as the main unit. *
No main control unit detected. Check system.	No main unit found. Check the system.	› Configure one device as the main unit. *
Auxiliary unit has error.	Fault on auxiliary unit	› Check the fault on the affected auxiliary unit *

* Only to be done by service technicians.

15 Transporting the unit



WARNING


Infection due to contaminated unit

- › Disinfect the unit before transport.
- › Close all media connections.



Wear protective equipment to avoid any risk of infection (e.g. liquid-tight protective gloves, protective goggles, face mask).

1. Prior to disassembly, clean and disinfect the suction unit and the unit via aspiration of a suitable disinfectant approved by the manufacturer.
2. Disinfect a defective unit using a suitable surface disinfection agent.
3. Seal all connections with sealing caps.
4. Pack the unit securely in preparation for transport.



- 1 Suction connection sealing cap (order number 7130-100-19)
- 2 Waste water connection sealing cap (order number 9000-412-98)
- 3 Sealing cap (order number 900041298)



Appendix

16 Menu structure

Setup wizard	Language				
	Operating mode				
	Exhaust air filter				
	Finished				
Start screen	Operating mode, vacuum				
	Performance, temperature				
	Control strategy options				
Settings	Access levels	Operator			
		Administrator			
			Service Technician [PIN]	XXXXXX	
	Device information	Device data	REF		
			SN		
			Firmware		
			Upgrade		
			Library version		
			PCB serial number		
			File		
				Device usage data	Number of starts
		System settings	Language	Deutsch, English, . . . ,	
	Date / time ²		Automatic		
			Date	DD MM YYYY	
			Time	HH MM	
			Time zone	UTC	
	Network ²		DHCP		
			IP address		
			Netmask		
		Gateway			
		MAC			
		Factory settings ¹	Clear message history		

Settings (continued)	Parameters ²	Cluster settings	Main unit
			Auxiliary unit 1
			Auxiliary unit 2
		Lag time ³	XX s
		Start / stop ³	Quick Start
			Eco Stop
		Exhaust air filter	Exhaust filter present
			No exhaust air filter
		Control strategy options ³	Adaptive operation
			Parallel operation
			Tyscor Tandem Flex
	Message history ²	List of alarms	Alarm information
	Maintenance	Filter maintenance	Maintenance completed
		Waste valve maintenance ²	Maintenance completed
		Exhaust air filter maintenance ²	Maintenance completed
	Service ²	Emergency mode	Reset
		Device components	Motor 0 (top)
			Motor 1 (bottom)
			Separation ¹

¹ only as Service Technician

² only as Administrator

³ only as Administrator and only on the main machine

17 Handover record

This document confirms that a qualified handover of the medical device has taken place and that appropriate instructions have been provided for it. This must be carried out by a qualified adviser for the medical device, who will instruct you in the proper handling and operation of the medical device.

Product name	Order number (REF)	Serial number (SN)

- ☐ Visual inspection of the packaging for any damage
- ☐ Unpacking the medical device and checking for damage
- ☐ Confirmation of the completeness of the delivery
- ☐ Instruction in the proper handling and operation of the medical device based on the operating instructions

Notes:

Name of person receiving instruction:

Signature:

Name and address of the qualified adviser for the medical device:

Date of handover:

Signature of the qualified adviser for the medical device:

--	--

18 Country representatives

Country	Address
UA  UA.TR.099	<p>Уповноважений представник в Україні: Приватне підприємство “Галіт” вул. 15 квітня, 6Є, с. Байківці, Тернопільський р-н, 47711, Україна тел.: 0800 502 998; +38 050 338 10 64 www.galit.te.ua; e-mail: office@galit.te.ua</p> <p>Виробник: Дюрр Дентал ЕсЕ Хьопфігхаймер Штрассе 17, Д-74321 Бітігхайм-Біссінген, Німеччина email: info@duerrdental.com</p>
CN	<p>备案人/生产企业：DÜRR DENTAL SE 德国迪珥齿科股份公司 住所/生产地址：Höfigheimer Str. 17, 74321 Bietigheim-Bissingen, Germany 联系方式： 电话：+ 49 7142 705-0 邮箱：info@duerrdental.com 网址：www.duerrdental.com</p> <p>代理人/售后服务单位：迪珥医疗器械（上海）有限公司 住所：上海市长宁区天山路 641 号 2 号楼 (20 幢) 303 室 联系方式： 电话：+ 86 21 6381 0270 传真：+ 86 21 6381 0290 邮箱：info@duerr.cn 网址：http://www.duerrdental.com</p>



Hersteller / Manufacturer:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
Fon: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

