

VistaSoft PatientBridge

EN Manual

Contents



Important information

1	About this document	2
---	-------------------------------	---



Installation

2	Requirements	3
3	Installing the software	3
4	Licensing	4
5	Uninstalling the software	4



Configuration

6	Creating a profile	5
6.1	Creating a profile	5
6.2	Configuring the profile	5
7	Delete profile	6



Operation

8	Launching PatientBridge	7
9	Importing patient data	7
10	Starting a data transmission directly	8

Important information

1 About this document

The following information is provided as a supplement to the VistaSoft manual (order-number: 2110100001).



<http://qr.duerredental.com/2110100001>



As a general rule, the VistaSoft handbook must also be consulted, as it contains important information relating to e.g. safety instructions, installation, configuration etc.

This manual forms part of the software. It corresponds to the version of the software and the technical standards valid at the time of installation.



In the event of non-compliance with the instructions and information contained in this manual, the manufacturer and the distributor will not offer any guarantee or accept any liability for the safe operation and the safe functioning of the unit and the software.

The German version of the instructions is the original manual. All other languages are translations of the original manual.



Installation

2 Requirements

PatientBridge has been developed to enable patient data from third-party software (e.g. accounting software) to be transferred to the imaging software.



It only supports third-party software that is based on either the framework *Win-Forms* or *WPF*.

3 Installing the software

Requirements:

- ✓ VistaSoft is installed
- ✓ Third-party software is installed
- › Open the installation package (VistaSoft_*.iso).
- › Open the ...\\Tools**PatientBridge** directory.
- › Start the **VistaSoft PatientBridge.msi** file.
- › Click **Next** to start the installation wizard.
- › Accept the Licence Agreement.
- › Specify the installation directory.
The path can be amended by clicking on **Change....**
- › Click **Finish** to complete the installation process.

4 Licensing

The software cannot be used without a valid license file.

Licensing is handled via VistaSoft and is described in the VistaSoft manual (order no.: 2110100001) in the section *Licensing*.

5 Uninstalling the software

The software is uninstalled via the control panel of the operating system:

- › Select the program in the *Control Panel*.
- › Click *Uninstall*.
- › Confirm uninstallation in the dialog window.

Configuration

VistaSoft PatientBridge Config opens automatically the first time PatientBridge is started up. To change the configuration, *PatientBridge Config* must be started manually:

- › Start configuration via **Start > All Programs > VistaSoft PatientBridge > VistaSoft PatientBridge Config**.

6 Creating a profile

A separate profile must be created for every third-party software from which data is to be transferred.

Once it has been created, each profile needs to be configured (see "6.2 Configuring the profile").



It only supports third-party software that is based on either the framework *Win-Forms* or *WPF*.

6.1 Creating a profile

- › Click **Create Profile...**


The *PatientBridge Profile Creation* window opens.

- › Enter a name under **Profile Name**.

- › Click **Select**.

The mouse cursor changes to .

- › Click the title bar of the third-party software.



The name of the third-party software appears in the **Window** line.

If the third-party software is not supported, the error message *Window is not supported!* is displayed. In this case the third-party software cannot be used with PatientBridge.

- › Click **Create** to save the profile.

6.2 Configuring the profile

Before data can be copied from the third-party software, the individual fields (e.g. patient name, date of birth) must be assigned.

- › Select the corresponding profile under **Select profile:**.

All fields that are yet to be mapped to fields from the third-party software are marked with




- › Select an entry (e.g. patient ID).

- › Under **Settings**, click **Select Patient ID**.

The mouse cursor changes to .

- › In the third-party software, click the field containing the patient ID.

This field must contain a value. Empty fields cannot be imported.




The *Property Selection* window opens.



If the corresponding field cannot be selected, it is either protected (content cannot be selected and copied), or the third-party software is not supported (see "2 Requirements").

- › Select the entry containing the *patient ID* from the list.

- › Check that the correct value has been applied in the configuration.



The field has been mapped successfully.

- › Proceed in the same way for the other fields.

i With some fields you will need to configure the corresponding notation or presentation of the information.

For example, under "Name" – is this stated as "Surname, Name" in a combined field, or are there separate fields for each?

For example, under "Date of birth": Is the date format in the third-party software DD.MM.YYYY or MM.DD.YYYY?

- › Click **Save** to save the profile.

7 Delete profile

- › Select the corresponding profile under **Select profile:**.
- › Click **Delete Profile...**.
- › In the dialog window, confirm with **Delete**.
The selected profile and all corresponding settings are deleted.



8 Launching PatientBridge

- › Start PatientBridge via *Start > All Programs > VistaSoft > VistaSoft PatientBridge* .

or

Double click the icon on the desktop.



9 Importing patient data

Before data can be transferred from the third-party software to VistaSoft, the following requirements must be met:

- ✓ PatientBridge is open
- ✓ VistaSoft is open
- ✓ Third-party software is open
- › Close any ongoing X-ray image acquisitions or tasks in VistaSoft.
 - Log out any patient who is logged in.
- › Bring PatientBridge to the foreground.
- › Select the required profile.



The last used profile is automatically activated in PatientBridge.

- › Click  to open the source window in the third-party software.
- › Open the patient in the third-party software for which you wish to transfer data.
- › Click .

The patient data is transferred to VistaSoft and the patient is automatically logged in.

10 Starting a data transmission directly

In the standard configuration, PatientBridge needs to be started and **Send** needs to be clicked for every transmission. Here, the PatientBridge window is permanently in the foreground. A link can be created as a way to get around this. Once this link is clicked, PatientBridge is called up and the data is transmitted once.

- ✓ Profile created in PatientBridge (see "6 Creating a profile")
- › Create a link to the file *PatientBridge.exe* – e.g. right-click on > Send to > Desktop (Create link).
- › Rename the link accordingly so that it is not confused with the original file.
- › Right-click the link and select *Properties*.
- › Under *Destination*: add */start* after the displayed path.
Examples:
"C:\Users\[user name]\Desktop\PatientBridge.exe" /start
C:\Users\[user name]\Desktop\PatientBridge.exe /start
- › Apply the changes and close the Properties window.

Whenever the link is activated, PatientBridge is launched in the background and the patient data is imported automatically.



As PatientBridge runs in the background, no error messages or similar are displayed.



Hersteller / Manufacturer:

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
Fon: +49 7142 705-0
www.duerrdental.com
info@duerrdental.com

