


# VistaSoft Cloud View



EN      Manual



# Contents



## Important information

<b>1</b>	<b>About this document</b>	<b>2</b>
1.1	Explanation of symbols used	2
1.2	Copyright information	2
<b>2</b>	<b>Safety</b>	<b>2</b>
2.1	Intended purpose	2
2.2	Improper use	2
2.3	General safety information	2
2.4	Protection from threats from the Internet	2



## Product description

<b>3</b>	<b>System requirements</b>	<b>3</b>
----------	----------------------------	----------



## Usage

<b>4</b>	<b>Notes on use</b>	<b>4</b>
<b>5</b>	<b>Starting the software</b>	<b>6</b>
<b>6</b>	<b>Open image</b>	<b>6</b>
<b>7</b>	<b>Display image information</b>	<b>7</b>
<b>8</b>	<b>Navigating</b>	<b>8</b>
8.1	Navigating in CBCT images (3D)	8
<b>9</b>	<b>Edit image temporarily</b>	<b>11</b>
9.1	3D	11
9.2	Display	18
9.3	Image corrections	18
9.4	Histogram	18
9.5	Measurements	19
9.6	Drawings	20
9.7	Annotations	20
<b>10</b>	<b>Close image</b>	<b>22</b>



## Troubleshooting

<b>11</b>	<b>Tips for operators and service technicians</b>	<b>23</b>
-----------	---	-----------



## Appendix

<b>12</b>	<b>Additional licensing terms</b>	<b>24</b>
12.1	Apache License 2.0	24
12.2	MIT License	27
12.3	DCMTK License	28
12.4	GLEW License	28
12.5	Mesa 3D License	29
12.6	BSD licence	30
12.7	FreeBSD licence	31

## Important information

### 1 About this document

This manual forms part of the software. It corresponds to the version of the software and the technical standards valid at the time of installation.



In the event of non-compliance with the instructions and information contained in this manual, the manufacturer and the distributor will not offer any guarantee or accept any liability for the safe operation and the safe functioning of the unit and the software.

The German version of the instructions is the original manual. All other languages are translations of the original manual.

#### 1.1 Explanation of symbols used

These symbols are used within this document and in the software:



Note, e.g. specific instructions regarding efficient and cost-effective use of the unit.



Manufacturer

#### 1.2 Copyright information

All circuits, processes, names, software programs and units mentioned in this document are protected by copyright.

The Installation and Operating Instructions must not be copied or reprinted, neither in full nor in part, without written authorisation from the copyright owner.

## 2 Safety

### 2.1 Intended purpose

The software offers only functions for displaying digital or digitalized 2D/3D X-ray images and videos. The software is exclusively for information and education purposes as well as for the use in regards to data exchange intended.

### 2.2 Improper use

The software must not be used for diagnostic purposes.

Information obtained with the help of this software must not be used for the purposes of any treatment or planning of a treatment.

The software does not have any specific filter for diagnosis or functions for managing patients. Images on which filters are used cannot be stored.

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damages resulting from improper usage. The user bears the sole risk.

### 2.3 General safety information

- › When using the software, always comply with all guidelines, laws and other rules and regulations that are applicable at the site of operation.
- › Do not modify the software.
- › Observe software manual.

### 2.4 Protection from threats from the Internet

The web application is downloaded from the Internet. The computer must therefore be protected from threats from the Internet.

- › Use antivirus software and update it regularly. Look for evidence of possible virus infection and, if applicable, check with the antivirus software and remove the virus.
- › Use, properly configure and regularly update your firewall.
- › Update the computers operating system regularly.



## Product description

VistaSoft Cloud View is a web application for looking at x-ray or video images.

The software is downloaded from the Internet and executed locally in the browser. The images remain in the local storage.

The images can be opened and modified with editing tools but not saved. The selection of editing tools depends on the opened image format.

The following image formats can be displayed:

- DICOM (\*.dcm)
- PNG (\*.png)
- JPEG (\*.jpg)
- TIFF (\*.tif)




The web application has not been tested or certified for clinical use. It is not approved as a medical device and must not be used for diagnostic purposes.

### 3 System requirements

Browser:	Google Chrome Microsoft Edge Opera (Each in the current PC version)
Interface:	Internet connection

## 4 Notes on use



- |   |                  |  |
|---|------------------|--|
| 1 | Patient data bar | Information on the patient and image if these are contained in the image file  |
| 2 | Light table      | Area in which the image is displayed   |
| 3 | Menu bar         | Various commands for the displayed image   |
| 4 | Toolbox          | Tools for temporary image processing<br>The tools on offer differ depending on the image type.   |
| 5 | Help             | Opens the software help in a separate browser tab  |
| 6 | Status bar       | The status bar displays notes on operation, e.g. on the use of the tools in the toolbox.<br>A new message pulses blue-grey for around 3 seconds. |

The menu bar contains various commands for the image displayed on the light table, depending on the image type:



Undo or repeat action



Hide/show annotations and measurements



Reset all MPR and TSA planes to their standard orientations (3D only)



Hide/show all MPR/TSA lines and the panoramic curve (3D only)



## 5 Starting the software

- › Open web application in the browser:



- › <https://view.vistasoft.com>

## 6 Open image

*Open image with drag & drop:*

- › Click and hold down the mouse button and drag the image file into the browser window (drag & drop).

- › Release the mouse button.

*Open image from the local file:*

- › Click on *Select images ....*
- › Select an image file and click on *Open*.



## 7 Display image information

The image information can be displayed for each image.

Unable to edit image information.


- › Click  in the patient data bar.

The image information included in the image is displayed in a flyout.

- › To change the display of the dental notation system (FDI, UNS, Palmer), select the desired dental notation system in the selection list.


## 8 Navigating

The following navigation aids can generally be used:

- Mouse wheel: zoom in/zoom out image section
- Right mouse button:  Move image section

### 8.1 Navigating in CBCT images (3D)

In digital volume tomography (CBCT), three-dimensional volume data is reconstructed from a large number of individual two-dimensional X-ray images (layer images). Various views are available in the software (see "Views").



*Fig. 1: Example: MPR view*

Various functions can be used to navigate through the slice planes/layers in the CBCT images. With the aid of these functions the slice planes can be moved or rotated in such a way that the important regions can be made visible for the examination.

In the illustration below all of the navigation functions have been made visible.













Fig. 2: Example image showing all of the navigation functions

	Jump to this point	Double-click in any view (except the volume view) to jump to this point in the volume. The slice planes in all views are moved accordingly.
	Move individual slice plane	Press and hold the left mouse button on a slice plane to drag it horizontally/vertically.
	Move two slice planes together	Press and hold the left mouse button on the intersection between two slice planes to move them together at the same time.
	Rotate slice plane	Press and hold the left mouse button on the dotted edge of a slice plane to rotate it.
	Tilt slice plane	Press and hold the left mouse button in the edge region of the view to tilt the view.
	Scroll layers	Press and hold the left mouse button in a free area in the view to scroll through the layers.

An orientation button  is displayed at the bottom right of each view showing the current viewing angle onto the patient.

The view can be changed in the volume view. The slice planes are displayed in the volume view, but they cannot be changed here.

	Move image section	Press and hold the right mouse button to move the image section upwards/downwards or to the left/right.
	Rotate view	Press and hold the left mouse button to rotate the view.
Zoom image section in/out		Rotate the mouse wheel to increase or reduce the image section. The zoom tools of the toolbox ( <i>Display</i> ) are not active in the volume view.

The orientation head always displays the current viewing angle onto the patient. In the volume view the viewing angle can be selected directly via the orientation head:

---



Orientation head front

Move the mouse cursor onto the orientation head. The views onto the orientation head are displayed and can be clicked on to select them.



Orientation head rear



Orientation head left



Orientation head right



Orientation head top



Orientation head bottom

## 9 Edit image temporarily

The toolbox contains tools for temporarily changing the presentation of the images on the light table.

The tools are divided into groups and can be expanded or collapsed with ►.

Selection of tools depends on which image is open on the light table.

Image changes can be reversed with ↶ and restored with ↷.



The changes to the image are not saved.

All changes are deleted when the image is closed.

### 9.1 3D

#### Views


In digital volume tomography (CBCT), three-dimensional volume data is reconstructed from a large number of individual two-dimensional X-ray images (layer images). Different views are available for evaluation of the data. These views can be selected under **View**.

##### Panoramic view



In the panoramic view a panoramic view is also displayed in addition to the slices. This panoramic view is calculated on the basis of the panoramic curve along the arch of the jaw. In addition, a slice of a transversal layer image (TSA slice) is displayed, which is calculated orthogonally to the panoramic curve.

Every CBCT image is initially opened in the panoramic view.




- 1 Panoramic view (calculated along the arch of the jaw)
- 2 Volume view
- 3 Sagittal MPR slice (view from left/right)
- 4 TSA slice (orthogonal to the panoramic curve)
- 5 Axial MPR slice (view from the top/bottom)

## TSA view



With the TSA view, several transversal layer images (TSA) are displayed in addition to the panoramic view.

The TSA slices are numbered (bottom left). The numbering allows adjacent TSA slices to be recognised, which makes drawing the nerve channel easier.




- 1 Panoramic view (calculated along the arch of the jaw)
- 2 Volume view
- 3 TSA slices (orthogonal to the panoramic curve), the number of slices can be chosen via the layout
- 4 Axial MPR slice (view from the top/bottom)

If the mouse cursor is positioned above one of the TSA slices, this TSA slice is highlighted in colour in the other views.

## MPR view









With the multiplanar reconstruction (MPR), any two-dimensional layer and slice images can be displayed.



- 1 Axial MPR slice (view from the top/bottom)
- 2 Volume view
- 3 Coronal MPR slice (view from the front/rear)
- 4 Sagittal MPR slice (view from left/right)

## Panoramic

	Standard panoramic curve, lower jaw	The standard curve form is automatically positioned on an axial layer that extends through the lower jaw
	Standard panoramic curve, upper jaw	The standard curve form is automatically positioned on an axial layer that extends through the upper jaw
	Standard panoramic curve on current slice plane	The standard curve form is manually positioned on the current axial slice. The centre point of the arch of the jaw is manually positioned by the user.
	User-defined panoramic curve on current slice plane	A free curve form is manually positioned on the current axial slice. The individual points of the panoramic curve are manually positioned one after the other.
	Abort creation of user-defined panoramic curve	During creation of the user-defined panoramic curve the points that have already been positioned are deleted again.
	Panoramic thickness	The layer thickness of the panoramic curve can be freely adjusted with a slider. With a lower thickness the panoramic view is sharper but noisier. With increasing thickness the noise is reduced, but the panoramic view appears less sharp. Here, with increasing thickness the panoramic view approaches the character of a panoramic X-ray image (2D).

The panoramic curve (except for the anatomically adjusted panoramic curve) can be moved and edited once it has been positioned.

The panoramic curve is also shown in the annotations.

### ***Positioning the panoramic curve automatically:***

- › Click  (maxilla) or  (mandible).


The panoramic curve is positioned. The panoramic view is calculated.

### ***Position standard panoramic curve on current slice plane:***

- › Adjusts the slice plane in the axial MPR slice.
- › Click .
- › In the axial MPR slice click on the centre point of the arch of the jaw.

The standard panoramic curve is positioned. The panoramic view is calculated.

### ***Creating a user-defined panoramic curve:***

- › Adjusts the slice plane in the axial MPR slice.
- › Click .
- › In the axial MPR slice left-click on the edge of the arch of the jaw.
- › Position the further points of the arch of the jaw with the left mouse button.
- › Place the end point of the arch of the jaw with the right mouse button.

The panoramic view is calculated.

### ***Moving a panoramic curve:***

- › Press and hold the left mouse button on the line of the panoramic curve (not on a point) and move the curve.





The entire panoramic curve is moved. The active layer is moved at the same time in parallel.

### ***Editing a panoramic curve:***


- › To move a point of the panoramic curve, press and hold the mouse button on the point and drag it.
- › To add a point, press the Ctrl key and left-click on the panoramic curve between the points.
- › To delete a point, press the Ctrl key and left-click on the point on the panoramic curve.
- › To scale the panoramic curve, press the Shift key, left-click on the panoramic curve and hold the button pressed while increasing/reducing the panoramic curve.



**TSA**

	Select TSA layout	Flyout for selection of the layouts or the number of transversal slices that are to be displayed (maximum of 6 rows x 10 columns)
	TSA distance	With the aid of a slider the distance between the transversal slices can be changed (0.2 mm - 10.0 mm).
	TSA length	With the aid of a slider the length of the transversal slices can be changed (2 mm - 200 mm).
	TSA width	With the aid of a slider the width of the transversal slices can be changed (0.1 mm - 10.0 mm).


**Select TSA layout:**

- › Click .
- The flyout with the matrix for selection of the TSA views opens.
- › Move the mouse cursor over the matrix.
- › When the required TSA views are selected (marked blue) click on the matrix.


**Result:**

The TSA slices are displayed and shown in the other views.


**Moving TSA slices:**

- › In the panoramic view or the axial slice, use the left mouse button  to move the TSA slices.
- All TSA slices are moved together along the panoramic curve.

**Tilting TSA slices:**



- › In the panoramic view, use the left mouse button  to tilt the TSA slices.
- The TSA slices can be tilted by max.  $\pm 75^\circ$ .
- TSA slices are tilted together in parallel. The position of the TSA slices is also displayed in the volume view, but cannot be changed in this view.


**MPR**

	Thickness of MPR slices	The thickness of the MPR slices can be changed with the slider. The lower the thickness, the sharper but the noisier the image. The higher the thickness, the lower the noise but the less sharp the image.
---	-------------------------	--

The adjustment always applies to all slices (sagittal, coronal, axial) of an MPR view.


## Nerve canal

	Draw nerve canal in image	The nerve channel can be drawn in all two-dimensional views (panoramic view and slices). The individual anchor points of the nerve channel are manually positioned one after the other. Once the nerve channel has been drawn, it is visible as a semitransparent "tube" in all views, including the volume view. It is saved as an annotation.
	Delete nerve channel	Deletes a nerve channel that has already been created. To do this, the nerve channel must be selected in a view or in the annotations.

 In the list of annotations, the drawn nerve channel can be called up again. When it is called up, the slices and the panoramic view are reset to the representation in which the nerve channel was drawn.

The nerve channel can be drawn in all two-dimensional views (slices and panoramic view) and is displayed in all views (including the volume view).

### *Drawing the nerve canal manually:*

- › Click .
- › Left-click with the mouse to position the anchor points of the nerve channel.  
You can also hold the shift key and use your mouse wheel to scroll through the layers.
- › Set the last anchor point with the right mouse button.

### *Edit nerve channel:*

- › To move an anchor point of the nerve channel, click on the anchor point, hold the mouse button and move the point.
- › To add an anchor point, press the Ctrl key and left-click on the nerve channel between the anchor points.
- › To delete an anchor point, press the Ctrl key and left-click on the point of the nerve channel.

### *Change display of nerve channel:*

- › To change the colour and thickness, select the nerve channel in the annotations and click .




### *Delete nerve channel:*

- › Select the nerve channel in a two-dimensional view or in the annotations.
- › Click .

## Volume

There are different display options in the volume view:

X-ray	Simulated X-ray projection image
Maximum value	Simulated X-ray projection image, the brightest point is shown for every image point in the viewing direction
ISO layer	View of the volume combining two opaque surfaces: soft tissue (brown) and bone structures (white)
3D	3D rendering
User-defined	3D rendering with simulated illumination. Skin tissue and muscle tissue are displayed in addition.

Slice		Display of the slice planes in the volume view
	Show slice plane edges only	Only the edges of the slice planes are displayed (may be masked by the volume display).
	Show slice planes with semitransparent surface	The slice planes are displayed with colour-filled surfaces (may be masked by the volume display or other slice planes).
	Hide slice planes	The slice planes are not visible.
Axial		Region above or below the slice plane is hidden.
Sagittal		Region to the left or right of the slice plane is hidden.
No slice		Everything is visible, all regions are displayed.
Reverse		Hidden regions are inverted (top/bottom or left/right).

In the **3D** view different tissue types are shown with different colours. The visibility of the tissue types can be adjusted with the relevant slider.

Tissue types in the **3D** view:

- Bone
- Dentine
- Enamel
- Metal

Tissue types in the **User-defined** view:

- Skin
- Muscle
- Bone
- Dentine
- Enamel
- Metal

The colours for the tissue types can each be chosen via a flyout in the **User-defined** view.

## 9.2 Display

These tools allow you to alter the image view. The zoom tools change the size that the image is displayed on the light table. The settings are not saved.

The settings for rotating and inverting an image are not saved.

With CBCT images (3D) the zoom tools can only be used on the MPR slices (axial, sagittal and coronal). The tools for image rotation and mirroring are not available for CBCT images (for information on how to rotate a slice, see "8.1 Navigating in CBCT images (3D)").

### Zoom



This matches the size of the image to the size of the window.



This magnifies the image so that one screen pixel equals one image pixel.



The current zoom factor is displayed. A defined zoom factor can be selected. (only 2D)



The zoom can be adjusted to any required setting via the slider. You can also use the mouse wheel to adjust the zoom. (only 2D)

### Rotate image



Rotate the image anti-clockwise by 90°.



Rotate the image clockwise by 90°.

### Mirroring



Horizontally invert the image



Vertically invert the image

## 9.3 Image corrections

The image brightness, contrast and gamma can be corrected. As you adjust the sliders, you can see the changes directly in the image itself.



Adjust the image brightness



Adjust the image contrast



Adjust the image gamma value  
< 1: Bright areas become darker

> 1: Dark areas become brighter

### Mouse pointer icon



Brightness and contrast can also be changed directly in the image using the left mouse button. For this purpose, go to **Toolbox > Image corrections** and enable the option **Mouse mode**. When the left mouse button is held down, making horizontal motions over the image changes the brightness, while vertical motions adjust the contrast. The sliders in the toolbox move along accordingly.

This option can be enabled separately for 2D X-ray/video images and for 3D X-ray images. It is preset for 2D X-ray/video images.


## 9.4 Histogram

The histogram shows the brightness distribution of the image.

Histogram windowing allows you to alter the image's grey scale values. The **black limit** shows the percentage of grey scale values changing to black, and the **white limit** shows the percentage of grey scale values changing to white. The limits can be altered to suit your requirements. The settings you have applied will be immediately apparent in the X-ray image on the light table.

The changes this makes to the image can emphasize structures that are otherwise too small to be distinguished by the human eye so that they can be identified.

### Editing a histogram:

- › Use the two sliders to adjust the black and white limits.
- › To move the entire selected range between the black and white limits click  in the middle of the region and move it with the mouse button pressed.

## 9.5 Measurements



In order to estimate lengths and angles in 2D X-ray images, the image needs to be calibrated with the aid of a reference image object. Despite calibration, this estimation does not show the accuracy of a measurement and must not be used as a measured result. The accuracy depends greatly on the projection distortion of the object on the image detector area.

In video images, lengths and angles cannot be estimated to scale. Panoramic X-ray images are not suitable for measuring due to their projection technique.

On CBCT images (3D) the calibration is calculated geometrically. Measurements are only possible in the slices (axial, sagittal, coronal and TSA), but not in the volume view or panoramic view.

The accuracy of the measurement results depends largely on the accuracy of the metadata generated by the software that exported the image.

The following measurements can be performed:



Simple line (start and end point)



Polyline (start, intermediate point and end point)



Angle (between two straight lines)

The image needs to be calibrated for the measurements.



Calibration  
Calibrate the image using this length for the selected line

Rulers can be displayed on calibrated images. Rulers are not available for uncalibrated images or for the measurement unit 'pixels' (px).




Show rulers



Hide rulers

**Perform calibration with a reference object:**


- › Click .
- › In the image, click on the start point with the left mouse button.

- › Click with the left mouse button on the end point.
- › Under **Tools** and **Calibration** enter the actual length of the measured length (e.g. the diameter of a steel sphere).




The values are copied to the image.


**Measuring a simple line:**

- › Click .
- › In the image, click on the start point with the left mouse button.
- › Click with the left mouse button on the end point.  
The measured result is displayed directly on the measurement line.


**Measuring a polyline:**

- › Click .
- › In the image, click on the start point with the left mouse button.
- › Use the left mouse button to add as many intermediate points as you want.
- › Click with the right mouse button on the end point.  
The measured result for the overall distance is displayed directly on the measurement line.

**Measuring an angle:**

- › Click .
- › In the image, click on the start point of the first straight line with the left mouse button.
- › Click with the left mouse button on the end point of the first straight line.
- › Click on the start point of the second straight line with the left mouse button.
- › Click with the left mouse button on the end point of the second straight line.  
The angle between the two straight lines is displayed.


**Editing a measurement:**

- › In the image, click on the measurement with the left mouse button.  
This activates the measurement and the anchor points are displayed.
- › To move the anchor point, click on the anchor point with the left mouse button , hold the mouse button and move it.

**Change the display of the measurement:**

- › Click the measurement to activate it.
- › Change the colour or line width.

**Delete measurement:**

- › Click the measurement to activate it.
- › Click .

## 9.6 Drawings

Drawings can be inserted into the image. All drawings are saved as annotations.

The following drawing tools are available:



Create free-hand drawing



Create a line



Create an arrow




Create a rectangle




Create an ellipse


### Create free-hand drawing:

- › Click .
- › In the image, click on the start point of the free-hand drawing using the left mouse button.
- › Keep left mouse button down and create a freehand drawing.
- › Release the mouse button to finish the free-hand drawing.  
The mouse cursor retains the function of the free-hand drawing.
- › To create a second free-hand drawing, press the left mouse button and keep it pressed while you draw.
- › To finish the free-hand drawing, press the right mouse button.


### Create a line:

- › Click .
- › In the image, click on the start point of the line using the left mouse button.
- › Click on the end point of the line using the left mouse button.

### Create an arrow:


- › Click .
- › In the image, click on the start point of the arrow using the left mouse button.
- › Click on the endpoint of the arrow using the left mouse button.

### Create a rectangle:


- › Click .
- › Click in the image with the left mouse button to place the first corner of the rectangle.

- › Move the mouse cursor.  
A preview of the rectangle is displayed.
- › Click with the left mouse button to define the size of the rectangle.

### Create an ellipse:

- › Click .
- › In the image, click on the start point of the ellipse using the left mouse button.
- › Move the mouse cursor.  
A preview of the ellipse is displayed.
- › Click with the left mouse button to define the size of the ellipse.

### Shift drawing:


- › In the image, click on the drawing using the left mouse button.  
This activates the drawing and the anchor points are displayed.
- › Keep the mouse cursor  pressed and move the drawing.
- › Release the mouse button.

### Edit drawing:


- › In the image, click on the drawing using the left mouse button.  
This activates the drawing and the anchor points are displayed.




The anchor points of pre-hand drawings cannot be edited.

- › To move the anchor point, click on the anchor point with the left mouse button , hold the mouse button and move it.

### Change the display of the drawing:

- › To change the colour and line thickness, select the drawing in *Annotations* and click .

### Delete drawing:

- › Mark the drawing in *Annotations* in order to activate it.
- › Click .

## 9.7 Annotations

All drawings, measurements and nerve channels of an image are saved as annotations. The anno-

tations are displayed in the list and can be called up again.

When an annotation is called up on a 3D image, the views/slices return to the view in which the annotations were drawn.

Every measurement is indicated in the list with the type, colour and numerical value.

### Type of annotation



Freehand annotation



Line drawing



Arrow drawing



Rectangle drawing



Ellipse drawing



Simple line measurement



Polyline measurement



Angle measurement



Nerve channel (3D only)



Panoramic curve




Implant

### Display annotation in the image:

- › Click on an annotation in the list.


If the annotation cannot be shown in the current view, the annotation cannot be selected.

### Change the display of an annotation:


- › Click  next to an annotation in the list.  
The flyout opens.
- › Select the colour.
- › For lines: select the line thickness.
- › For nerve channels: adjust the diameter of the line.

- › Click outside the flyout to close it.


### Edit annotation:

- › Click on an annotation in the list or an annotation in the image.  
The anchor points of the annotation are visible.
- › Move the mouse cursor over one of the anchor points of the annotation.  
The mouse cursor changes to .
- › Press and hold the left mouse button and drag the anchor point to the required position.
- › Release the mouse button.

### Shift annotation:

- › Click on an annotation in the list or an annotation in the image.  
The anchor points of the annotation are visible.
- › Click on the annotation between the anchor points.  
The mouse cursor changes to .
- › Press and hold the left mouse button and drag the annotation to the required position.
- › Release the mouse button.

### Delete annotation:

- › Click .
- This deletes the annotation. The deletion cannot be undone.




The annotations are not stored in the image. All annotations are deleted when the image is closed.

## 10 Close image



All changes are deleted when the image is closed.

- › In the browser, click on  *Back* to open another image.  
or  
Close browser window.



# ? Troubleshooting

## 11 Tips for operators and service technicians

Error	Possible cause	Remedy
<b>Presentation of the software surface is not correct</b>	Software not optimised for this browser	› Use another browser, see "3 System requirements".
<b>Presentation of 3D images too slow</b>	Graphic card does not meet the requirements	› Activate hardware acceleration in the browser. › Use a faster graphic card.

## Appendix

### 12 Additional licensing terms

Parts of the product are based on third-party image libraries with a large number of available licences that are not covered by the EULA.

Those that require a reproduction of the license text in the distribution are listed below.

- "12.1 Apache License 2.0"

- "12.2 MIT License"

- "12.3 DCMTK License"

- "12.4 GLEW License"

- "12.5 Mesa 3D License"

- "12.6 BSD licence"

- "12.7 FreeBSD licence"

Parts of the product are based on third-party software licensed under the GNU Lesser General Public License v2.1 (LGPL-2.1):

- FTGL library

- LibUsbDotNet library

- Qt framework

Parts of the product are based on the work of the Independent JPEG Group.

Some free software licences request but do not require confirmation of software use. We would like to express our gratitude to the following people and projects for their contributions to this product:

- zlib compression library (Jean-loup Gailly, Mark Adler and Team)

- libpng graphics library (Glenn Randers-Pehrson and Team)

- LibTIFF library (Glenn Randers-Pehrson and Team)

- SQLite database engine (D. Richard Hipp and Team)

### 12.1 Apache License 2.0

This license applies to the following redistributed libraries and files:

- Castle Windsor
  - Castle.Windsor.dll
- Castle.Core.dll
- Castle.Facilities.Logging.dll
- Castle.Services.Logging.Log4netIntegration.dll
- log4net: Copyright © 2004-2013 Apache Software Foundation
  - log4net.dll
- log4cxx: Copyright © 2012 Apache Software Foundation
  - log4cxx.dll
- Lucene.Net: Copyright © 2013 Apache Software Foundation
  - Lucene.Net.dll
  - Lucene.Net.Contrib.Analyzers.dll
  - Lucene.Net.Contrib.Highlighter.dll
  - Lucene.Net.Contrib.Memory.dll
- ASP.NET Core: Copyright © .NET Foundation and Contributors. All rights reserved.
  - Various DLLs
- Vaadin Components: Copyright © 2021 Vaadin Ltd. All rights reserved.
  - Included in dd-components.js
- Material Design Web Components: Copyright © 2018 Google Inc.
  - Included in dd-components.js

Apache License

Version 2.0, January 2004

<http://www.apache.org/licenses/>

TERMS AND CONDITIONS FOR USE,  
REPRODUCTION, AND DISTRIBUTION

1. Definitions.

"License" shall mean the terms and conditions for use, reproduction, and distribution as defined by Sections 1 through 9 of this document.  
"Licensor" shall mean the copyright owner or entity authorized by the copyright owner that is granting the License.

"Legal Entity" shall mean the union of the acting entity and all other entities that control, are controlled by, or are under common control with that entity. For the purposes of this definition, "control" means

(i) the power, direct or indirect, to cause the direction or management of such entity, whether by contract or otherwise, or (ii) ownership of fifty percent (50%) or more of the outstanding shares, or (iii) beneficial ownership of such entity.

"You" (or "Your") shall mean an individual or Legal Entity exercising permissions granted by this License.

"Source" form shall mean the preferred form for making modifications, including but not limited to software source code, documentation source, and configuration files.

"Object" form shall mean any form resulting from mechanical transformation or translation of a Source form, including but not limited to compiled object code, generated documentation, and conversions to other media types.

"Work" shall mean the work of authorship, whether in Source or Object form, made available under the License, as indicated by a copyright notice that is included in or attached to the work (an example is provided in the Appendix below).

"Derivative Works" shall mean any work, whether in Source or Object form, that is based on (or derived from) the Work and for which the editorial revisions, annotations, elaborations, or other modifications represent, as a whole, an original work of authorship. For the purposes of this License, Derivative Works shall not include works that remain separable from, or merely link (or bind by name) to the interfaces of, the Work and Derivative Works thereof.

"Contribution" shall mean any work of authorship, including the original version of the Work and any modifications or additions to that Work or Derivative Works thereof, that is intentionally submitted to Licensor for inclusion in the Work

by the copyright owner or by an individual or Legal Entity authorized to submit on behalf of the copyright owner. For the purposes of this definition, "submitted" means any form of electronic, verbal, or written communication sent to the Licensor or its representatives, including but not limited to communication on electronic mailing lists, source code control systems, and issue tracking systems that are managed by, or on behalf of, the Licensor for the purpose of discussing and improving the Work, but excluding communication that is conspicuously marked or otherwise designated in writing by the copyright owner as "Not a Contribution."

"Contributor" shall mean Licensor and any individual or Legal Entity on behalf of whom a Contribution has been received by Licensor and subsequently incorporated within the Work.

2. Grant of Copyright License. Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable copyright license to reproduce, prepare Derivative Works of, publicly display, publicly perform, sublicense, and distribute the Work and such Derivative Works in Source or Object form.

3. Grant of Patent License. Subject to the terms and conditions of this License, each Contributor hereby grants to You a perpetual, worldwide, non-exclusive, no-charge, royalty-free, irrevocable (except as stated in this section) patent license to make, have made, use, offer to sell, sell, import, and otherwise transfer the Work, where such license applies only to those patent claims licensable by such Contributor that are necessarily infringed by their Contribution(s) alone or by combination of their

Contribution(s) with the Work to which such Contribution(s) was submitted. If You institute patent litigation against any entity (including a cross-claim or counterclaim in a lawsuit) alleging that the Work or a Contribution incorporated within the Work constitutes direct or contributory patent infringement, then any patent licenses granted to You under this License for that Work shall terminate as of the date such litigation is filed.

4. Redistribution. You may reproduce and distribute copies of the Work or Derivative Works thereof in any medium, with or without modifications, and in Source or Object form, provided that You meet the following conditions:

(a) You must give any other recipients of the Work or Derivative Works a copy of this License; and  
(b) You must cause any modified files to carry prominent notices stating that You changed the files; and

(c) You must retain, in the Source form of any Derivative Works that You distribute, all copyright, patent, trademark, and attribution notices from the Source form of the Work, excluding those notices that do not pertain to any part of the Derivative Works; and

(d) If the Work includes a "NOTICE" text file as part of its distribution, then any Derivative Works that You distribute must include a readable copy of the attribution notices contained within such NOTICE file, excluding those notices that do not pertain to any part of the Derivative Works, in at least one of the following places: within a NOTICE text file distributed as part of the Derivative Works; within the Source form or documentation, if provided along with the Derivative Works; or, within a display generated by the Derivative Works, if and wherever

such third-party notices normally appear. The contents of the NOTICE file are for informational purposes only and do not modify the License. You may add Your own attribution notices within Derivative Works that You distribute, alongside or as an addendum to the NOTICE text from the Work, provided that such additional attribution notices cannot be construed as modifying the License.

You may add Your own copyright statement to Your modifications and may provide additional or different license terms and conditions for use, reproduction, or distribution of Your modifications, or for any such Derivative Works as a whole, provided Your use, reproduction, and distribution of the Work otherwise complies with the conditions stated in this License.

5. Submission of Contributions. Unless You explicitly state otherwise, any Contribution intentionally submitted for inclusion in the Work by You to the Licensor shall be under the terms and conditions of this License, without any additional terms or conditions. Notwithstanding the above, nothing herein shall supersede or modify the terms of any separate license agreement you may have executed with Licensor regarding such Contributions.

6. Trademarks. This License does not grant permission to use the trade names, trademarks, service marks, or product names of the Licensor, except as required for reasonable and customary use in describing the origin of the Work and reproducing the content of the NOTICE file.

7. Disclaimer of Warranty. Unless required by applicable law or agreed to in writing, Licensor provides the Work (and each Contributor provides its Contributions) on an "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND,

either express or implied, including, without limitation, any warranties or conditions of TITLE, NON-INFRINGEMENT, MERCHANTABILITY, or FITNESS FOR A PARTICULAR PURPOSE. You are solely responsible for determining the appropriateness of using or redistributing the Work and assume any risks associated with Your exercise of permissions under this License.

8. Limitation of Liability. In no event and under no legal theory, whether in tort (including negligence), contract, or otherwise, unless required by applicable law (such as deliberate and grossly negligent acts) or agreed to in writing, shall any Contributor be liable to You for damages, including any direct, indirect, special, incidental, or consequential damages of any character arising as a result of this License or out of the use or inability to use the Work (including but not limited to damages for loss of goodwill, work stoppage, computer failure or malfunction, or any and all other commercial damages or losses), even if such Contributor has been advised of the possibility of such damages.

9. Accepting Warranty or Additional Liability. While redistributing the Work or Derivative Works thereof, You may choose to offer, and charge a fee for, acceptance of support, warranty, indemnity, or other liability obligations and/or rights consistent with this License. However, in accepting such obligations, You may act only on Your own behalf and on Your sole responsibility, not on behalf of any other Contributor, and only if You agree to indemnify, defend, and hold each Contributor harmless for any liability incurred by, or claims asserted against, such Contributor by reason of your accepting any such warranty or additional liability.

END OF TERMS AND CONDITIONS

## 12.2 MIT License

This license applies to the following redistributed libraries and files:

- Json.Net:
  - Newtonsoft.Json.dll
- JsonCpp:
  - Various DLLs
- .NET: Copyright © .NET Foundation and Contributors. All rights reserved.
  - Various DLLs
- .NET Framework: Copyright © .NET Foundation and Contributors. All rights reserved.
  - Various DLLs
- Electron: Copyright © Electron contributors, Copyright © 2013-2020 GitHub Inc.
  - Various DLLs and files
- date-fns: Copyright © 2020 Sasha Koss and Lesha Koss
  - Included in dd-components.js
- nodeclrhst: Copyright © 2019 Daniel Martin
  - Various DLLs and files
- ONNX Runtime: Copyright © Microsoft Corporation
  - onnxruntime.dll
  - onnxruntime\_providers\_shared.dll
- DirectML: Copyright © Microsoft Corporation. All rights reserved.
  - DirectML.dll
- zxcvbn-ts: Copyright © 2012-2016 Dan Wheeler and Dropbox, Inc.

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF

MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NON-INFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

## 12.3 DCMTK License

This license applies to the redistributed DCMTK libraries, including the following files:

- dcmtk.dll

Copyright (C) 1994-2019, OFFIS e.V.  
All rights reserved.

This software and supporting documentation were developed by  
OFFIS e.V. R&D Division Health  
Escherweg 2

26121 Oldenburg, Germany

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.

- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

- Neither the name of OFFIS nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR

CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

## 12.4 GLEW License

This license applies to the redistributed GLEW library (OpenGL Extension Wrangler Library) file:

- glew32.dll

The OpenGL Extension Wrangler Library

Copyright (c) 2002-2008, Milan

Ikits

Copyright (c) 2002-2008, Marcelo E. Magallon

Copyright (c) 2002, Lev Povalahev

All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- \* Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.

- \* Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

- \* The name of the author may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE

LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

## 12.5 Mesa 3D License

This license applies to the redistributed Mesa 3-D Graphics Library file:

– opengl32.dll

Copyright (c) 1999–2007 Brian Paul.  
All Rights Reserved.

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL BRIAN PAUL BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

## 12.6 BSD licence

This license applies to the following redistributed libraries and files:

- FluentNHibernate: Copyright © 2008-2018, James Gregory and contributors. All rights reserved.
  - FluentNHibernate.dll
- lit-element: Copyright © 2021 The Polymer Authors. All rights reserved.
  - Included in dd-components.js
- Polymer: Copyright © 2017 The Polymer Authors. All rights reserved.
  - Included in dd-components.js
- OpenCV:
 

Copyright © 2000-2020, Intel Corporation, all rights reserved.

Copyright © 2009-2011, Willow Garage Inc., all rights reserved.

Copyright © 2009-2016, NVIDIA Corporation, all rights reserved.

Copyright © 2010-2013, Advanced Micro Devices, Inc., all rights reserved.

Copyright © 2015-2016, OpenCV Foundation, all rights reserved.

Copyright © 2015-2016, Itseez Inc., all rights reserved.

Copyright © 2019-2020, Xperience AI, all rights reserved.

Third party copyrights are property of their respective owners.

  - All distributed DLLs having a name starting with "opencv\_"
- OpenMesh: Copyright © 2001-2019, RWTH-Aachen University, Department of Computer Graphics and Multimedia. All rights reserved.
  - Included in ImagingLib.dll

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
3. Neither the name of the copyright holder nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.



## 12.7 FreeBSD licence

This license applies to the following redistributed libraries and files:

- OpenJPEG

This license applies to the redistributed OpenJPEG binary file: openjp2.dll

The copyright in this software is being made available under the 2-clauses BSD License, included below. This software may be subject to other third party and contributor rights, including patent rights, and no such rights are granted under this license.

Copyright © 2002-2014, Universite catholique de Louvain (UCL), Belgium

Copyright © 2002-2014, Professor Benoit Macq

Copyright © 2003-2014, Antonin Descampe

Copyright © 2003-2009, Francois-Olivier Devaux

Copyright © 2005, Herve Drolon, Freemage Team

Copyright © 2002-2003, Yannick Verschuere

Copyright © 2001-2003, David Janssens

Copyright © 2011-2012, Centre National d'Etudes Spatiales (CNES), France

Copyright © 2012, CS Systemes d'Information, France

All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE

COPYRIGHT HOLDER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.







**Hersteller / Manufacturer:**

DÜRR DENTAL SE  
Höpfigheimer Str. 17  
74321 Bietigheim-Bissingen  
Germany  
Fon: +49 7142 705-0  
[www.duerrdental.com](http://www.duerrdental.com)  
[info@duerrdental.com](mailto:info@duerrdental.com)

